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NEOXONE® - Steril PDS (Polidioksanon) Cerrahi İplik

Tanım

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NEOXONE® - Sterile PDS (Polydioxanone) Surgical Suture

Description

NEOXONE® PDS, is a sterile, synthetic and absorbable monofilament surgical suture material obtained by polymerization of the p-dioxanone monomer. It is dyed in FDA approved (D&C Violet No.2) to make it easily distinguishable. NEOXONE® PDS has a reliable and predictable absorption profile, as well as having very high initial tensile and knot-pull strength. Characteristics like its monofilament structure, smooth surface, highly uniform thread radius, minimal suture memory, high tensile strength, softness and elasticity provides excellent knot placement and tie-down, high knot safety and perfect handling, as well as smooth passage through tissues with minimal resistance, thereby inflicting minimal trauma on tissues. NEOXONE® PDS meets all requirements of the European Pharmacopoeia (EP - Pharm. Eur.) and the United States Pharmacopoeia (USP).

Indications

NEOXONE® PDS is indicated for use in general soft tissue approximation and/or ligation, where middle-long term wound support is required. NEOXONE® PDS is suitable for use in applications like orthopedics, plastics surgery, gynecology, gastrointestinal tract surgery and pediatric cardiovascular surgery, as well as being suitable for use on patients with longer than normal wound healing times because of lower blood supply levels and/or other underlying causes, like obese or old patients or patients with cancer or diabetes. It's not intended for use in cardiovascular or neurological surgery and thus should not be used in those procedures.

Mode of Action

NEOXONE® PDS elicits a minimal inflammatory reaction in tissues where it is used, which is a typical endogenous reaction against a foreign substance. As time passes, the suture material is encapsulated in fibrous connective tissue. NEOXONE® PDS is broken down by hydrolysis and is completely absorbed by metabolism, without causing any permanent changes to the wound site. Absorption begins

with loss of mass which is followed by a loss of tensile strength. After the day of its application, NEOXONE® PDS retains approximately 70% - 80% of its initial tensile strength on 14th day post implantation and 50% - 60% of its initial tensile strength by the 28th day. It provides effective wound support for up to 8-9 weeks. Absorption is minimal in the first 90 days and in normal perfusion of tissues, NEOXONE® PDS's absorption is completed between 180 - 210 days. It is highly resistant to bacterial colonization due to its monofilament and homogenous structure. Also, the reaction it elicits in tissues is minimal because of its physical and chemical properties.

Contraindications

NEOXONE® PDS is contraindicated for use in cases where permanent wound support by the suture material is required (like cardiovascular surgery and neurosurgery). NEOXONE® PDS should not be used in patients with known sensitivities and allergies to its components.

Warnings

NEOXONE® PDS should not be resterilized. Opened or damaged packages should not be used! Unused open packs should be properly discarded. NEOXONE® PDS should be kept in a cool and dry place (between 10°C-30°C), away from humidity and/or direct sunlight. Prolonged exposure of suture packs to extreme temperatures should be avoided. Should only be used by or under the order / supervision of a qualified medical doctor.

Notes / Precautions

While handling this suture material, extreme care should be taken to avoid damaging the suture thread by crushing, kinking and/or crimping it with surgical instruments such as forceps, needle holders, tweezers etc. Users should possess adequate knowledge on appropriate suture selection and suturing techniques, as the risk of experiencing wound dehiscence and other similar negative outcomes will vary with factors like tissues that the sutures are applied, suitable

suture selection and user's qualification. NEOXONE® PDS should be used in accordance with the generally accepted surgical practices and techniques of its application domain. In certain cases, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Before applying this suture in tissues with weak blood supply, it should be considered that late absorption and suture extrusion may occur. Sutures that remain longer than 7 days on skin can cause local irritation and should be removed or snipped off. In cases where quick healing of the wound is desired, the sutures should be removed after the required period has elapsed, without waiting for the absorption to be completed.

Adverse effects

Reported adverse effects include slight local irritation on the skin, slight tissue inflammation and wound dehiscence. Like all other suture materials, prolonged contact with factors such as urine and bile can lead to lithiasis.

How supplied

NEOXONE® PDS is available in sizes ranging from USP 8/0 (0.4 metric) to USP 2 (5 metric), either undyed or colored in Violet, with needles and in shrink wrapped boxes of a dozen of sterile, double layer unit packs. NEOXONE® PDS is sterilized with Ethylene Oxide. The needles are made of high quality 300 or 400 series stainless steel. To enable smooth passage through tissues, the needles are coated with medical grade silicone.

CE 2292

Symbols and signs used in labelling

STERILE EO Sterilized with Ethylene Oxide

LOT Lot no

Manufacturing date (MM YYYY)

Expiry date (MM YYYY)

For single use only

Caution! Check for specific warnings or precautions

Consult instructions for use

Keep dry and away from moisture

Keep away from sunlight

Do not use damaged, open or torn single packs

Keep/store between 10°C-30°C

Keep/store below 65% relative humidity

Do not resterilize

Manufacturer

Suture in loop form (double thread)

Needle thickness in microns (for microsurgery needles)

Heavy (thick) needle

MANUFACTURER: SETPA TIBBİ GEREÇLER İTH. İHR. SAN. VE TİC. LTD. ŞTİ.
11454 SOK. NO:13 KAT:1 YENİŞEHİR - İZMİR / TÜRKİY
TEL: +90 232 458 90 51 FAX: +90 232 469 38 77 - info@setpa.com.tr
www.setpa.com.tr