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NEOPLENE® Ultra - Steril POLİPROPİLEN Cerrahi İplik

Tanım

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NEOPLENE® Ultra - Sterile POLYPROPYLENE Surgical Suture

Description

NEOPLENE® Ultra is a sterile, synthetic and monofilament non-absorbable surgical suture material made of polypropylene. It is dyed in blue or black to make it easily distinguishable. NEOPLENE® Ultra's monofilament structure, ideal tensile strength, smooth and slippery surface, minimal pack memory, extra smoothness and controlled linear elongation provides excellent knot placement and tie-down, high knot safety and perfect handling, as well as providing smooth passage through tissues with minimum resistance, thereby inflicting minimal trauma on tissues. NEOPLENE® Ultra meets all requirements of the European Pharmacopoeia (EP - Pharm. Eur.) and the United States Pharmacopoeia (USP).

Indications

NEOPLENE® Ultra is indicated for use in general soft tissue approximation and/or ligation, where long term wound support - and therefore use of non-absorbable sutures - is required. NEOPLENE® Ultra is particularly suitable for use in ophthalmic surgery, orthopedics, general skin closures, plastics surgery, cardiovascular surgery and neurological surgery.

Mode of Action

NEOPLENE Ultra elicits a minimal inflammatory reaction in tissues where it is used, which is a typical endogenous reaction against a foreign substance. NEOPLENE® Ultra cannot be broken down by enzymes or hydrolysis in-vivo, it's biologically inert against degradatory and weakening effects of enzymes. As time passes, the suture material is encapsulated in fibrous connective tissue. Due to its monofilament and homogenous structure, NEOPLENE® Ultra is resistant to bacterial colonization so it is the ideal suture for use in contaminated or infected wounds. Having minimal interaction with tissues, it is useful as a pull-out suture (eg. subcuticular closure)

Contraindications

NEOPLENE® Ultra is contraindicated for use in patients who are known to be highly sensitive or allergic to polypropylene. Sutures that remain longer than 7 days on skin can cause local irritation and should be removed or snipped off. NEOPLENE® Ultra should not be used in patients with known sensitivities and allergies to its components.

Warnings

NEOPLENE® Ultra should not be resterilized. Opened or damaged packages should not be used! Unused open packs should be properly discarded. NEOPLENE® Ultra should be kept in a cool and dry place (between 10°C-30°C), away from humidity and/or direct sunlight. Prolonged exposure of suture packs to extreme temperatures should be avoided. Should only be used by or under the order / supervision of a qualified medical doctor.

Notes / Precautions

While handling this suture material, extreme care should be taken to avoid damaging the suture thread by crushing, kinking and/or crimping it with surgical instruments such as forceps, needle holders, tweezers etc. Users should possess adequate knowledge on appropriate suture selection and suturing techniques, as the risk of experiencing wound dehiscence and other similar negative outcomes will vary with factors like site of application, suitable suture selection and user's qualification. NEOPLENE® Ultra should be used in accordance with the generally accepted surgical practices and techniques of its application domain. In certain cases, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon.

Adverse effects

Reported adverse effects include slight local irritation on the skin, slight tissue inflammation and wound dehiscence. Like all other suture materials, prolonged contact with factors like urine and bile can lead to lithiasis. It may cause mild pain, edema and erythema at the wound site.

How supplied

NEOPLENE® Ultra is available in sizes ranging from USP 10/0 (0.2 metric) to USP 2 (5 metric), colored in Blue or Black, in needed form and in shrink wrapped boxes of a dozen of sterile, double layer unit packs. NEOPLENE® Ultra is sterilized with ethylene oxide. The needles are made of high quality 300 or 400 series stainless steel. To enable smooth passage through tissues, the needles are coated with medical grade silicone.

Symbols and signs used in labelling

- Sterilized with Ethylene Oxide
- Lot no
- Manufacturing date (MM YYYY)
- Expiry date (MM YYYY)
- For single use only
- Caution! Check for specific warnings or precautions
- Consult instructions for use
- Keep dry and away from moisture
- Keep away from sunlight
- Do not use damaged, open or torn single packs
- Keep/store between 10°C-30°C
- Keep/store below 65% relative humidity
- Do not resterilize
- Manufacturer
- Suture in loop form (double thread)
- Needle thickness in microns (for microsurgery needles)
- Heavy (thick) needle

MANUFACTURER: SETPA TIBBI GEREÇLER İTH. İHR. SAN. VE TİC. LTD. ŞTİ.
11454 SOK. NO:13 KAT:1 YENİŞEHİR - İZMİR / TÜRKİY
TEL: +90 232 458 90 51 FAX: +90 232 469 38 77 - info@setpa.com.tr
www.setpa.com.tr