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NEOSILK® Ultra - Steril İPEK Cerrahi İplik

Tanım

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NEOSILK® Ultra - Sterile SILK Surgical Suture

Description

NEOSILK® Ultra is a sterile, braided (twisted in 8/0 and finer sizes) and non-absorbable surgical suture material that is composed of an organic protein called fibroin, derived from 100% natural silk fibers obtained from the cocoons of silk worms by a special purification process called degumming. It is dyed in black to make it easily distinguishable and has successfully passed the dye permanence tests performed by TSE. NEOSILK® Ultra is coated with silicone (except for 8/0 and finer sizes) which provides non-capillarity and enables the suture to easily and reliably slide without damaging itself while tying knots. It also provides smooth passage through tissues with minimal resistance, thereby inflicting minimal trauma on tissues. Furthermore, NEOSILK® Ultra's ideal tensile strength, smooth and slippery surface obtained by silicone coating, minimal pack memory and extra smoothness provides excellent knot placement and tie-down, high knot safety and perfect handling. NEOSILK® Ultra meets all requirements of the European Pharmacopoeia (EP - Pharm. Eur.) and the United States Pharmacopoeia (USP).

Indications

NEOSILK® Ultra is indicated for use in all patients for general soft tissue approximation and/or ligation, where long term wound support - and therefore use of non-absorbable sutures - is required, including ophthalmic surgery. NEOSILK® Ultra is also suitable for use in applications like general skin closures, gastrointestinal tract surgery, plastics surgery, cardiovascular surgery and neurological surgery.

Mode of Action

NEOSILK® Ultra elicits a minimal inflammatory reaction in tissues where it is used, which is a typical endogenous reaction against a foreign substance. As time passes, the suture material is encapsulated in fibrous connective

tissue. While being classified as a non-absorbable suture, silk suture gradually loses its tensile strength, hence its wound support, due to the gradual in-vivo breaking up of the organic protein fibers that it is made of. It is resistant to bacterial colonization due to its smooth surface and non-capillarity obtained by silicone coating.

Contraindications

NEOSILK® Ultra is contraindicated for use in patients who are known to be highly sensitive or allergic to silk. Since the initial tensile strength of silk is gradually lost in-vivo, it should not be used in cases where permanent wound support by the suture material is required. It is also contraindicated in infected wounds and in urinary and biliary tract surgery. NEOSILK® Ultra should not be used in patients with known sensitivities and allergies to its components.

Warnings

NEOSILK® Ultra should not be resterilized. Opened or damaged packages should not be used! Unused open packs should be properly discarded. NEOSILK® Ultra should be kept in a cool and dry place (between 10°C-30°C), away from humidity and/or direct sunlight. Prolonged exposure of suture packs to extreme temperatures should be avoided. Should only be used by or under the order / supervision of a qualified medical doctor.

Notes / Precautions

While handling this suture material, extreme care should be taken to avoid damaging the suture thread by crushing, kinking and/or crimping it with surgical instruments such as forceps, needle holders, tweezers etc. Users should possess adequate knowledge on appropriate suture selection and suturing techniques, as the risk of experiencing wound dehiscence and other similar negative outcomes will vary with factors like site of application, suitable suture selection and user's qualification.

NEOSILK® Ultra should be used in accordance with the generally accepted surgical practices and techniques of its application domain. In certain cases, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Sutures that remain longer than 7 days on skin can cause local irritation and should be removed or snipped off.

Adverse effects

Reported adverse effects include slight local irritation on the skin, slight tissue inflammation, enhanced bacterial activity and wound dehiscence. Like all other suture materials, prolonged contact with factors like urine and bile can lead to lithiasis. It may cause mild pain, edema and erythema at the wound site.

How supplied

NEOSILK® Ultra is available in sizes ranging from USP 8/0 (0.4 metric) to USP 2 (5 metric), colored in Black (also in Blue for USP 8/0), with or without needles and in shrink wrapped boxes of a dozen of sterile, double layer unit packs. NEOSILK® Ultra is sterilized with ethylene oxide. The needles are made of high quality 300 or 400 series stainless steel. To enable smooth passage through tissues, the needles are coated with medical grade silicone.

CE 2292

Symbols and signs used in labelling

Sterilized with Ethylene Oxide

Lot no

Manufacture date (MM YYYY)

Expiry date (MM YYYY)

For single use only

Caution! Check for specific warnings or precautions

Consult instructions for use

Keep dry and away from moisture

Keep away from sunlight

Do not use damaged, open or torn single packs

Keep/store between 10°C-30°C

Keep/store below 65% relative humidity

Do not resterilize

Manufacturer

Suture in loop form (double thread)

Needle thickness in microns (for microsurgery needles)

Heavy (thick) needle

MANUFACTURER: SETPA TIBBİ GEREÇLER İTH. İHR. SAN. VE TİC. LTD. ŞTİ.
11454 SOK. NO:13 KAT:1 YENİŞEHİR - İZMİR / TÜRKİYE
TEL: +90 232 458 90 51 FAX: +90 232 469 38 77 - info@setpa.com.tr
www.setpa.com.tr