

# Гемоблок® Haemoblock®

## LOCAL HEMOSTATIC IN LIQUID FORM



### FAST AND RELIABLE HEMOSTASIS


#### MEDICAL DEVICE

Haemoblock® is a local topical  
regenerating hemostatic agent in  
liquid form with an antiseptic effect


### A WIDE RANGE OF APPLICATIONS

Haemoblock® is used to stop all types of capillary and parenchymal bleeding that occur during surgical operations and injuries (including operations on parenchymal organs and bleeding from cancellous bone), for the treatment of wounds, burns, trophic ulcers, as well as for hemostasis and visualization during endoscopic manipulations, maxillofacial and ophthalmic surgery.

Before applications Haemoblock®



Upon the application of Haemoblock®



Prosthetics of the thoracoabdominal aorta

### APPLICATION METHOD

Haemoblock® — stops capillary, diffuse  
and parenchymal bleedings  
of any localization and at any stage  
of surgical operation

Swabs, generously soaked in  
Haemoblock® onto pre-dried  
wound surface



Irrigating the wound surface  
without follow-up rinsing off  
at any operation stage



Endoscopically (via the  
instrument channel)




HEMOSTASIS ACHIEVED  
WITHIN 2 MINUTES

**NOT APPROPRIATE  
FOR ARTERIAL / JET HAEMORRHAGES!**


### MECHANISM OF ACTION

Haemoblock® forms a complex with  
blood proteins

Haemostasis within 2 minutes



The resulting clot is lysed naturally



Microphotograph of a clot formed Haemoblock®

### ADVANTAGES

Reliable hemostasis behind 2 minutes  
Healing with primary suture tension without  
suppuration

Neutral pH  
Does not form scab, non-cauterizing,  
no recurring bleeding

Antiseptic effect due to  
Silver

Fast tissue regeneration  
Due to healing without the scab formation  
and the reduction of necrosis zone

Liquid form  
Use in endoscopy, laparoscopy, thoracoscopy  
(introduced into cavities through a laparoscopic  
port)

Transparent  
Full visualization of the wound surface



### DESCRIPTION

**MEDICAL DEVICE**  
**Haemostatic solution**  
**HAEMOBLOCK®**



RELEASE FORM:  
15 ml, 20 ml, 30 ml, 50 ml

**Composition:**  
Aqueous solution of the integral silver salt of  
Poly(Acrylic) acid, with integrated silver  
nanoparticles.

**Indications for use:**  
Diffuse, capillary and parenchymal bleedings;  
surgical operations, including parenchymal organs,  
injuries, wounds, burns.

REGISTRATION CERTIFICATE  
for a medical device No. FS2012/13587 dated November 16, 2021  
Haemostatic agent "Haemoblock" according to TU 9391-002-680873-2012



PRODUCTION:  
"PUL-SAR" Ltd, Russia

WhatsApp/Telegram  
+7 909 943 30 77  
info@pul-sar.com  
+7 906 062 06 18  
info@pul-sar.ru



[www.haemoblock.ru/en](http://www.haemoblock.ru/en)