

Actolind[®] w Solution


Cleaning, irrigation, moistening and care solution for wounds

THE OPTIMUM COMBINATION OF REFERENCE AND SMART MOLECULES FOR **WOUND CARE**


A **Synergistic** Approach in Wound Care with
Polyhexanide + Poloxamer

Your best choice for wound care

- Cleansing, moistening and decontamination of wounds
- Convenient and painless cleaning
- Removes scabs, necrotic tissues and biofilm layers from wounds effectively
- Does not inhibit granulation or epithelialization
- Suitable for antiseptic cleansing, softening and removal of encrusted dressing materials without causing a secondary trauma
- Highly effective against multidrug-resistant organisms (MRSA, VRE) and mycotic pathogens (*C. albicans*)
- Dermatologically tested and rated as “**excellent**”
- Biocompatible
- Well tolerated



Take care of the wounds with Actolind[®]



Actolind® w Solution

Actolind® w Solution is a cleansing, irrigation, decontamination, moistening and care solution used together with materials such as bandages, dressings, pads and wound fillers for acute, chronic and infected wounds as well as 1st and 2nd degree burns.

Actolind® w Solution cleans wounds mechanically. Surface active agents in its composition provide a convenient and painless cleaning needed before wound care by decreasing the surface tension. **Actolind® w Solution** contributes to the wound healing process by ensuring a highly positive cleansing and moisturizing effect on applied areas.

Actolind® w Solution is effective against multidrug-resistant (e. g. MRSA, VRE) and mycotic pathogens. It is also well tolerated by allergics and diabetics.

Chemical Composition:

Poloxamer, Polyhexamethylene biguanide (Polyhexanide, 0.1%) purified water

Indications:

- Infected or uninfected acute wounds (e.g. traumatic wounds that occur as the result of cuts, tears, excoriation, crushing, being bitten, being stung or surgical wounds)
- Infected or uninfected chronic wounds (diabetic ulcers, venous ulcers, decubitus (pressure sore), etc.)
- Infected or uninfected thermal or chemical wounds (1st and 2nd degree burns, chemical burns, wounds that caused by electrical current, radiation wounds etc.)
- Entry regions of urological catheters,
- Maintenance of peristomal skin, fistulas and abscesses, wound surface prior to dressing, suitable for antiseptic cleansing, softening and removal of dressing materials, bandages, dressings, pads and wound fillers that are dried, contaminated and adhered to the wound

Actolind® w Solution should be used for wound surface cleaning and irrigation before any treatment. **Actolind® w Solution** can be applied in combination with **Actolind® w Gel**. The wound should be cleaned with **Actolind® w Solution** before the application of **Actolind® w Gel** to remove any contaminants and deposits.

Actolind® w Solution can also be used for antiseptic cleansing and softening of encrusted dressing materials that are dried and adhered to the wound surface and difficult to remove.

Actolind® w Solution does not inhibit granulation or epithelialization. It has neither irritating nor pain giving properties. It causes no sensitization and can be used safely for a long time.

Usage:

Actolind® w Solution can be applied directly onto the wound as well as onto the bandage and gauzes.

The entire wound and its surroundings should be thoroughly cleaned with **Actolind® w Solution**. For optimal wound cleaning, it is recommended to wet the wound consisting of deep and intensively necrotic tissues, for at least 10-15 min. with the help of a fully wetted compress (application of wet dressing). The application can be repeated according to the wound conditions without any time limits. Do not dilute the solution! The cleaned wound can be covered with **Actolind® w Gel** and bandaged properly.

Thoroughly wet dressings and bandages adhered to the wound surface with **Actolind® w Solution**, wait 5-10 min. for softening and remove them carefully from wound. Wetting the dressing materials with **Actolind® w Solution** prevents adhering to the wound surface and facilitates dressing change significantly. In order to maintain a moist wound climate, the dressing can be changed daily, if necessary, several times a day. **Actolind® w Solution** is compatible with all modern wound dressings.

Possible infections can be avoided by regular cleaning catheter inlet ports, flexible connectors and its environment by **Actolind® w Solution**. Spray **Actolind® w Solution** on the area from a short distance and wait 1-5 min. The ingredients of **Actolind® w Solution** are particularly mucous membrane-friendly, cleanse catheter and PEG material without damaging it. As with wound care, encrusted sticking plasters can be removed easily after wetting with **Actolind® w Solution**.

Actolind® w Solution can be heated to body temperature prior to use.

Packaging			
Art. No	Unit	Packaging Type	Packaging Unit
06.6703.05	5 ml	Monodose Strip	100 x 5
06.6703.10	10 ml	Monodose Strip	50 x 5
06.6703.50	50 ml	Bottle with spray cap	200
06.6703.100	100 ml	Bottle with spray cap	100
06.6703.250	250 ml	Bottle with spray cap	40
06.6703.500	500 ml	Bottle	20
06.6703.1	1000 ml	Bottle	10

Do not use Actolind® w Solution in the following cases:

1. In hyaline cartilage and aseptic joint surgeries, in the middle and internal ear and in the eye (If **Actolind® w Solution** contacts the aseptic cartilage, wash immediately with Ringer solution or physiological saline solution)
2. In case of hypersensitive to any of the ingredients
3. On central nervous system or meninx
4. Along with anionic surfactants, cleansing soaps, creams, oils, enzymes, etc. (These substances should be removed fully from the wound surface before use)
5. Peritoneal lavage and washing

This product is a Class III Medical Device according to the Medical Device Directive 93/42/EEC



Actolind[®] w Gel


Protective moisturizing and care gel for acute,
chronic and infected wounds

THE OPTIMUM COMBINATION OF REFERANCE AND SMART MOLECULES FOR **WOUND CARE**


A **Synergistic** Approach in Wound Care with
Polyhexanide + Poloxamer

Your best choice for wound care

- Cleansing, moistening and decontamination of wounds
- Convenient and painless cleaning
- Removes scabs, necrotic tissues and biofilm layers from wounds effectively
- Does not inhibit granulation or epithelialization
- Suitable for antiseptic cleansing, softening and removal of encrusted dressing materials without causing a secondary trauma
- Highly effective against multidrug-resistant organisms (MRSA, VRE) and mycotic pathogens (*C. albicans*)
- Dermatologically tested and rated as “**excellent**”
- Biocompatible
- Well tolerated



Take care of the wounds with Actolind[®]



Actolind[®] w Gel

Actolind[®] w Gel is a product used for mechanical cleaning, decontamination and moistening of the chronic and infected wounds as well as 1st and 2nd degree burns together with medical devices such as bandage, dressing, pad and wound fillers.

Actolind[®] w Gel cleanses, moisturizes and disconnects the wound surface from outer environment by forming a protective film layer. It fills the space between dressing and wound surface and helps the healing process get shorter by keeping the wound stable.

Chemical Composition:

Poloxamer, Polyhexamethylene biguanide (Polyhexanide, 0.1%) water and auxiliary substances

Indications:

- Infected or uninfected acute wounds (e.g. traumatic wounds that occur as the result of cuts, tears, excoriation, crushing, being bitten, being stung or surgical wounds)
- Infected or uninfected chronic wounds (diabetic ulcers, venous ulcers, decubitus (pressure sore), etc.)
- Infected or uninfected thermal or chemical wounds (1st and 2nd degree burns, chemical burns, wounds that caused by electrical current, radiation wounds etc.)
- Entry regions of urological catheters
- Maintenance of peristomal skin, fistulas and abscesses, wound surface prior to dressing, suitable for antiseptic cleansing, softening and removal of dressing materials, i.e. bandages, dressings, pads and wound fillers that are dried, contaminated and adhered to the wound.

Actolind[®] w Gel can be applied in combination with **Actolind[®] w Solution**. The wound should be cleaned with **Actolind[®] w Solution** before the application of **Actolind[®] w Gel** to remove the contaminants and deposits. The wound bed should be covered with a plenty of gel prior to the application of dressings and bandages.

Actolind[®] w Gel does not inhibit granulation or epithelialization.

It has neither irritating nor pain giving properties. It causes no sensitization and can be used safely for a long time.

Actolind[®] w Gel can stay on the wound surface until the next dressing change. Thus, **Actolind[®] w Gel** keeps the wound surface constantly humid and minimizes the adhesion of the dressing on the wound surface by preventing of drying out.

Actolind[®] w Gel is effective against multidrug-resistant (MRSA, VRE) and mycotic pathogens. It is also well tolerated by allergics and diabetics. The protective film layer, formed by **Actolind[®] w Gel**, can be easily removed from the wound surface during the next dressing change. For this purpose, **Actolind[®] w Solution** is recommended. This application should be repeated on a regular basis to remove the dead tissues and scabs during the healing process.

Usage:

Actolind[®] w Gel can be applied directly onto the wound as well as onto the bandage and gauzes. Materials such as bandage, gauze, wiping clothes, etc. should not be used again. The opened product should be used on one single patient to avoid the risk of any contamination.

The entire wound and its surroundings should be thoroughly cleaned with **Actolind[®] w Solution** to minimize the risk of transmission of contaminants and microorganisms to the wound surface. For optimal wound cleaning, it is recommended to wet the wound, consisting of deep and intensively necrotic tissues, for at least 10-15 min. with the help of a fully wetted compress (application of wet dressing). Do not dilute the solution! The cleaned wound should be covered with **Actolind[®] w Gel** and bandaged properly.

Packaging			
Art. No	Unit	Packaging Type	Packaging Unit
06.6803.05	5 ml	Monodose Strip	100 x 5 strip
06.6803.10	10 ml	Monodose Strip	50 x 5 strip
06.6803.30	30 ml	Tube	240
06.6803.40	40 ml	Tube	200
06.6803.50	50 ml	Tube	200
06.6803.100	100 ml	Tube	200

Do not use Actolind[®] w Gel in the following cases:

1. In hyaline cartilage and aseptic joint surgeries, in the middle and internal ear and in the eye (If **Actolind[®] w Gel** contacts the aseptic cartilage, wash immediately with Ringer solution or physiological saline solution)
2. In case of hypersensitive to any of the ingredients
3. On central nervous system or meninx
4. Along with anionic surfactants, cleansing soaps, creams, oils, enzymes, etc. (These substances should be removed fully from the wound surface before use)

This product is a Class III Medical Device according to the Medical Device Directive 93/42/EEC

