HNK-NBL-S

Kompresörlü Nebulizatör

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Hünkar Ecza ve Medikal San. Tic. Ltd. Sti

2

KS/01/19.11.2018

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21

KS/01/19.11.2018

HNK-NBL-S Compressor Nebulizer User Manual

Hünkar Ecza ve Medikal San. Tic. Ltd. Şti

Index

1. Product Description	3
2. Product Specialities	3
3. Device Technical Specification	4
4. Device's Design and Introduction	5
5. Warnings and Measure of Safety	6
6. Device Installation and Usage Instructions	7
7. Cleaning and Disinfection	10
8. Maintenance and Clues	11
9. Device's Technical Datas	12
10. Keeping and Transporting	14
11. Troubleshooting	15
12. Warranty Conditions and After Sales Support	15
13. EMC (Electromagnetic Compatibility) Decalaration	16
14. Symbols and Descriptions	20
15. Device & Outside Box's Symbols Descriptions	21
16. European Conformity Mark	21
17. Manufacturer and Technical Service Information	21

This product is medical device. Before using the device please read user manual carefully and use it with in the direction of your doctor.

1. Product Description

Thanks for preferring **Nimo**® compressor nebulizer. **Nimo**® compressor nebulizer is designed for treatment of asthma, allergy and other respiration disease. HNK-NBL-S, converts electricity energy in to mechanical energy with it's compressor motor and this provides high pressure air outlet which converts the medicine in medication cup in to aerosol that respirable easily.

This device is suitable for patients from all different ages but it is not suitable for patients who is unconscious or cant breath himself. Babies, small children and adults who need special care should use it under the control of care-taker.

Nimo® compressor nebulizer gives the medication to patient silently, fastly and also without disturbing the patient. HNK-NBL-S is easy to use, light and ergonomic and is designed for use at home or at work.

2. Product Specialities

- It provides high density aerosol with it's powerful compressor motor.
- Low medication residual
- Light and portable compact design.
- Gives you freedom of action with it's long air tube.

3. Device Technical Specification


Product Name	Compressor Nebulizer
Model Number	HNK-NBL-S
Method of Operation	Compressor Electricity Motor
Nominal Operation Voltage	220 V AC 50Hz
Power Consumption	60 W
Nebulization Rate	>0.3 ml/min
MMAD	4.7±30%µm
Operating Pressure Range	0.6 – 1.1 bar (9 – 16 psi), (62 – 110 Kpa)
Compressor Pressure Range	2.4 – 3.4 bar (35 – 50 psi), (241 – 345 Kpa)
Liter Flow Range	4-6 liter/min
Medication Cup Capacity	Max 8 ml
Working Environment	Temperature: 5 °C – 40 °C Humidity: <90% RH Atmospheric pressure: 86-106 kPa
Storage and Transporting Conditions	Temperature: -25°C-70°C Humidity: 10-95% RH Atmospheric pressure: 70-106 kPa
Classification	Class II Equipment
Type	Type B Applied Parts
Ingress Protection	IP0X (Unguarded against dust or liquid)
Noise	<60 db (less than 1 mt)
Size	(width x length x height) 20 cm x 27 cm x 10 cm
Weight	1.4 kg

4

KS/01/19.11.2018

4. Device's Design and Introduction

a) Device's appearance and pieces

Piece Number	Piece Name	Piece Number	Piece Name
1	Device Main Unit	6	Medication Cup
2	Electricity Cable	7	Mouthpiece
3	Device ON/OFF Key	8	Air Tube
4	Air Filter Lid	9	Adult and Child Mask
5	Air Out Hole	10	Replacement Air Filters

b) Key Specialities

Key	Function
0	Device Off
1	Device On

5

KS/01/19.11.2018

5. Warnings and Measure of Safety

- ⚠ Please read user manually carefully before using the device.
- ⚠ The device is designed for patients from every age but babies, small children and adults who need special care should use it under the control of care-taker. Please keep away the device from children and babies. This device consists small pieces which can go down easily.
- ⚠ This device shouldn't use in the enviroment where inflammable gases, oxygen and aerosol sprays used in.
- ⚠ This device is used for breathing passage treatment under the control of doctor Please do not misuse it except for its purpose. If you feel yourself bad while using the device stop using and ask your doctor for help.
- ⚠ Patients who is unconscious and can't aspirate himself should't use this device.
- ⚠ Do not use in anaesthetic or ventilator breathing circuits.
- ⚠ Please follow your doctor's or respiration specialist advise for medication type and amount.
- ⚠ Please use original parts and accessories. Warranty is invalid if you use not original parts and user's personal reasons.
- ⚠ Before use the device be sure all the accessories are working Each user should use her/his own medication cup, mask and mouthpiece.
- ⚠ Please unplug after using the device and do not plug and unplug the device with wet hand.
- ⚠ Do not operate the device that fall in to water. Do not pour water or any liquids on to the device if happens unplug the device as soon as possible and dry the device.
- ⚠ Because of the drowning danger keep away air-tube and cable of the device from children.
- ⚠ Never close the air holes of the device or put the main unit the place where air holes can be closed.


6

KS/01/19.11.2018


6. Device Installation and Usage Instructions

⚠ Before first usage of the device the nebulizer and it's accosaries should be cleaned and disinfection. Check Chapter 7. "Cleaning Disinfection"

- 1) First open the whole package and take-out the device and accessories.
- 2) Install the air-tube to air out hole on the device as it is shown in the picture.

- 3) Install the other edge of air-tube to down side of medication cup as it is shown in the picture.


7

KS/01/19.11.2018


- 4) Turn the medication cup lid reverse of clockwise and pull up and take out it.
- 5) Put the right type and amount of medicine in to medication cup as your doctor advised to you.

⚠ While putting the medicine do not pass the maximum limits of medication cup. Maximum capacity of medication cup is 8 ml.

- 6) Put the medication cup on its place and turn clockwise the lid to close it.

- 7) As your usage choice install the mask or mouthpiece in to medication cup.


- 8) Be sure the air-filter is mounted on device and power key is "0" mode.

8

KS/01/19.11.2018

- 9) Plug the device in to socket which has suitable voltage and frequency.
- 10) Operate the device by pushing On/Off button of it and turn it "1" mode.
- 11) Be sure aerosol is coming from medication cup.
- 12) Wear the mask as it is shown in the picture and cover your mouth and nose. Keep breathing normally.

- 13) Put the mouthpiece between your lips and teeth as it is shown in the picture. Keep breathing normally.

⚠ Keep the medication cup upright otherwise nebulisation can be inefficient.

⚠ Be calm and relax in nebulisation process. Keep breathing normally.

- 14) Continue your inhalation as your doctor's advised duration or till medication finishes in medication cup.
- 15) When you finish inhalation push the On/ Off key and turn it in to "0" and unplug the nebulizer.

9

KS/01/19.11.2018

7. Cleaning and Disinfection

⚠ It can cause infection risk if you don't clean the medication cup regularly.

⚠ Clean the medication cup, medication cup lid, mouthpiece and masks after every usage. Disinfection should be done after the last usage of the day.

The special methods for cleaning and disinfection are given below:

Cleaning after every usage

Wash the medication cup, mouthpiece and mask in hot water which is not boiled and clean the pieces. (Less than 50°C)

⚠ Don't wash the device's main unit, air tube and air filter.

- 1) Take-out the device's Mouthpiece, mask, medication cup and air tube.
- 2) Empty the residual in medication cup.
- 3) Wash the medication cup, mask and mouthpiece in hot water.
- 4) Dry the washed parts with clean and soft cloth. Leave the parts on clean and dry cloth to get rid of humid and water on it.
- 5) Wipe the device's body and air tube with clean and humid cloth.
- 6) After cleaning all the parts be sure all the parts are dried.

⚠ Don't wash the parts in washing machine.

⚠ Don't dry parts in micro-wave.

Disinfection:

It is important disinfect the medication cup, medication cup lid, mouthpiece and mask after the last usage of the day.

- 1) Follow the after every usage cleaning steps to clean medication cup, medication cup lid, mouthpiece and mask.
- 2) Put the medication cup, medication cup lid, mouthpiece and mas in to boiled water and wait for 10 minutes.

10


KS/01/19.11.2018

8. Maintenance and Clues

- In normal usage conditions the air-filter should change every 60 days. The air filters which are blocked and dirtied too much can cause negative effects on nebulisation.

Changing Filter

- 1) First turn-off the device and unplug.
- 2) Take out the filter lid by pulling as it is shown in the picture. After lid is taken-out if filter stays in the device use tweezers to take-out.

- 3) Install the new filter in to lid.
- 4) Install the lid which filter is installed to its place.

11

KS/01/19.11.2018

- ⚠ Don't reuse the used filter by cleaning again.
- ⚠ Just use the filters which are produced by manufacturer otherwise the device can be harmed.
- ⚠ Never use device without filter.
 - Because of using frequently and disinfection processes of nebulizers and accessories plastic parts exposure fray. This can change the aerosol specialties which is coming from medication cup and effects the inhaling negatively. For this reason it is recommended to change medication cup, masks, air-tube and mouthpiece after 1 year.
 - Please protect the device, medication cup and accessories against hard impulse. Otherwise this can cause damage and device will become useless.
 - Clean the device and its parts before every usage.
 - Do not extract, fix or any change on device otherwise this can cause electric shock, fire or leakage.
 - Never use the device it is damaged or fall in to water.
- ⚠ If any problem occurs while using the device first check the troubleshooting. If problem continuous apply to the authorised service.

9. Device's Technical Datas

- The device is manufactured in compliance with EN 13544-1:2007+A1:2009 standards for airway treatment devices.
- The nebulisation system datas are given in the table in down side.


MMAD (Particle size)	4.7±30%µm
Nebulisation Rate	>0.3 ml/min
GSD (Geometric Standard Deviation)	4.77
Particle Size Under 5 µm	~%50

- The horizontal axis is the particle size value, the value is logarithmic distribution;

12

KS/01/19.11.2018

The left vertical axis is the cumulative percentage of the volume, corresponding to the rising trend of the curve.

- The small particle size (MMAD) for 50% cumulative particle mass indicates that the particle sizes at the aerosol from the nebulizer outlet can reach the lower airways.
- In the nebulizer, the particles in the range of 5 µm constitute approximately 50% of the aerosol output and more.
- GSD shows aerodynamic particle size distribution and is calculated as 4.77 for Nebulizer.

13

KS/01/19.11.2018

10. Keeping and Transporting

- The device can be kept with its package in well aired room humid is not more than %10-95 and the temperature between -25-70 °C.
- Under this storage condition explained up this device can be kept for 10 years.
- While carrying protect the package from getting wet, impacts and being crushed.

⚠ Keep the device away from children.

⚠ Keep the device away from direct sunlight and abrasive gases.

Package Contents:

- HNK-NBL-S Compressor Nebulizer.
- User Manual.
- Warranty.

Accessories inside the case

- Medication Cup (1 Pcs)
- Adult Mask (1 Pcs)
- Child Mask (1 Pcs)
- Mouthpiece (1 Pcs)
- Air Tube (1 Pcs)
- Spare Air-Filter (5 Pcs)

14

KS/01/19.11.2018

11. Troubleshooting

Problem	Solutions
Device's On/Off key is 1 position but doesn't work.	☞ Check the power cable is plugged or not. ☞ Check the device's body is overheated or not. If overheated too much unplug and wait for cooling.
Low aerosol or there is no aerosol.	☞ Check the medication cup is filled with water soluble and true medication. ☞ Check the medication cup is filled up with right size of medicine. ☞ Check the connection between device medication cup and tube is correct or not. ☞ Check the medication cup has right connection to device. ☞ Check air-tube has no blockage. ☞ Medication can be blockaged. Clean it as it is explained in chapter 7 "Cleaning and Disinfection" ☞ Check the medication cup is in vertical position or not. ☞ Check the air-tube is straight or not. ☞ The device can be clogged air filter, please check.
Device overheats too much	☞ Check the air channels on the device is clear and open. ☞ Don't cover the device with any cloths. ☞ In long inhaling process break for some time and then continue inhaling.

⚠ If these suggestions doesn't solve the problem please don't try to fix it and ask help from authorized service.

12. Warranty Conditions and After Sales Support

The device's warranty duration is 2 years after purchasing.

- ! If you follow the instructions in manual and charge and discharge it properly the durable period will be 10 years for this device.
- ! **Warning: The warranty will be invalid if the device is harmed, fixed or opened by user.**
- ! **Note: Please do not forget to take the warranty and sales receipts and keep them.**

15

KS/01/19.11.2018

13. EMC (Electromagnetic Compatibility) Decalaration

This equipment generates, uses and can radiate radio frequency (RF) energy. This equipment may cause electromagnetic interference if it is not installed and used in accordance with the manual.

This device has been tested for Medical Products in accordance with EN-60601-1-2 and has been determined to comply with acceptable limits. These limits indicate that the device provides an acceptable level of protection against electromagnetic interference (EMC) if the device is used as specified in the manual.

This device may be affected by portable and mobile RF communication devices. This equipment must not be stored together with other equipment.

See the following tables for more information about this device and EMC.

Guidance and Manufactures Declaration-Electromagnetic		
These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.		
Immunity Test	Compliance	Electromagnetic environment—guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufactures Declaration-Electromagnetic			
These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.			
Immunity Test	IEC60601 Test Level	Compliance Level	Electromagnetic environment—guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	± 6 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (UT de >=95 düşüş) 0,5 devir için %40 UT (UT de %60 düşüş) 5 devir için %70 UT (UT de %30 düşüş) 25 devir için <5 % UT (UT de >=95 düşüş) 5 saniye için	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power source
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the AC mains voltage prior to application of the test level.

Guidance and Manufactures Declaration-Electromagnetic			
These devices are intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that the device is used in such an environment.			
Immunity Test	IEC60601 Test Level	Compliance Level	Electromagnetic environment—guidance
Conducted RF IEC 61000-4-6	3 V rms 150 kHz ila 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter Recommended separation distance $d = \frac{1}{40} \sqrt{P}$ $d = 1,20 \sqrt{P}$ $d = 1,20 \sqrt{P}$ 80 MHz ila 800 MHz $d = 1,20 \sqrt{P}$ 800 MHz ila 2,5 GHz Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10 V rms 150 kHz ila 80 MHz 10 V/m 80 MHz ila 2,5 GHz	10 V 10 V/m	

- Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.
 Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
 a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.
 b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m

Recommended separation distance between portable and mobile RF communication equipment and this device				
These devices are intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	150 kHz ila 80 MHz, Outside ISM Bands	150 kHz ila 80 MHz, In ISM Bands	80 MHz ila 800 MHz	80 MHz ila 800 MHz
	$d = 1,16 \sqrt{P}$	$d = 1,20 \sqrt{P}$	$d = 0,4 \sqrt{P}$	$d = 7,64 \sqrt{P}$
0,01	0,12	0,12	0,12	0,23
0,1	0,37	0,38	0,38	0,73
1	1,16	1,20	1,20	2,30
10	3,67	3,79	3,79	7,27
100	11,60	12,00	12,00	23,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
Note 1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies..
Note 2. ISM (industrial, scientific, and medical) bands from 150 kHz to 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz.
Note 3. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

14. Symbols and Descriptions

	This device complies with the requirements of the Medical Devices Directive(93/42/EEC)
	Manufacturer
	Date of Manufacture
	Class II Equipment per IEC 60601-1
	Type B Equipment per IEC 60601-1
	Consult instructions for use
	LOT Number
	Model Number
	Serial Number
	Disposal of Electrical & Electronic Equipment (WEEE): Do not treat this product as household waste.

15. Device & Outside Box's Symbols Descriptions

	Warning/Caution/Note
	Keep Away From Sunlight
	Keep Away From Getting Wet
	Handle With Care
	Rohs Declaration
	Barcode Number

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The product conforms to European Medical Device Directive (93/42/EEC) and meets applicable health, safety and environmental requirements. If the mark is accompanied by a number, conformity is verified by the indicated notified body.

17. Manufacturer and Technical Service Information

Hünkar Ecza ve Medikal San. Tic. Ltd. Şti.
Akçeşme Mahallesi 2019 sk. No:5 Merkezefendi
Denizli /TÜRKİYE
Tel: +90 258 371 46 56
www.hunkarmedikal.com