

HALKALI ÖĞÜTÜCÜ

VIBRATORY DISC MILL



Orta-sert, sert, kırılğan ve zorlu malzemelerin çok hızlı, güvenli ve verimli bir şekilde öğütülmesi için kullanılır. Laboratuvarlarda ve endüstriyel tesislerde numune hazırlama adına oldukça kullanışlıdır. Uygulama Alanları: Alaşımlar, bazalt, çimento klinker, seramik, seramik tozu, kömür, kok kömürü, inşaat malzemeleri, feldspat, cam ve malzemeleri

ÜRÜN TEKNİK ÖZELLİKLERİ

Elektrik motoru: 1,5 HP, 1000 dev/vak, 380 v, 3 p, 50 Hz

50, 100 ve 250 cc kapasiteli 3 farklı öğütme seti seçeneği

Ürün Besleme Boyutu: < 15 mm*

Besleme Miktarı : 35 - 150 ml

Ürün Çıkış İnceliği: < 38 μ m*

Öğütme süresi 0'dan 999 saate, 999 dakikaya ya da 999 saniyeye kadar ayarlanabilir.

The Vibratory Disc Mill is used for the very rapid, safe and efficient grinding of medium-hard, hard, brittle and tough materials. It is a very useful for sample preparation in laboratories and industrial plants. Application Examples: Alloys, basalt, cement clinker, ceramics, chamotte, coal, coke, construction materials, feldspar, glass and various materials.

TECHNICAL DETAILS

Electrical Motor: 1,5 HP, 1000 dev/vak, 380 v, 3 p, 50Hz

Grinding Jar Sizes : 50cc, 100cc, 250cc

Feed Quantity: 35 - 150 ml

Feeding Size: < 15 mm*

Final Fineness: < 38 μ m*

Dijital setting of grinding time from 0 up to 999 h or 999 min. or 999 seconds

