


# TEST IT GEL




**BLOOD GROUPING**  
Gel Cards & Manual Systems


## GEL CARDS




**TEST-IT Gel**  
AB (x4)  
REF TG101  
CE 1434




**TEST-IT Gel**  
ABO/2D & Reverse  
REF TG102  
CE 1434




**TEST-IT Gel**  
ABO/CDE (DVI+) & Reverse  
REF TG103  
CE 1434




**TEST-IT Gel**  
ABO/CDE (DVI-) & Reverse  
REF TG104  
CE 1434




**TEST-IT Gel**  
ABO/Rh(C,E) DVI+ & Reverse  
REF TG105  
CE 1434




**TEST-IT Gel**  
ABO/Rh(C,E) DVI- & Reverse  
REF TG106  
CE 1434




**TEST-IT Gel**  
ABO/2D + Kell & Reverse  
REF TG107  
CE 1434




**TEST-IT Gel**  
Confirm ABO/DVI+  
REF TG108  
CE 1434




**TEST-IT Gel**  
Confirm ABO/DVI-  
REF TG109  
CE 1434




**TEST-IT Gel**  
Newborn  
REF TG110  
CE 1434




**TEST-IT Gel**  
DAT Screen  
REF TG111  
CE 1434




**TEST-IT Gel**  
Cross Match  
REF TG112  
CE 1434




**TEST-IT Gel**  
Type / Screen  
REF TG113  
CE 1434




**TEST-IT Gel**  
Coombs  
REF TG114  
CE 1434




**TEST-IT Gel**  
Neutral  
REF TG115  
CE




**TEST-IT Gel**  
Neutral / Coombs  
REF TG116  
CE 1434




**TEST-IT Gel**  
Rh Phenotype  
REF TG117  
CE 1434




**TEST-IT Gel**  
Rh Phenotype + Kell  
REF TG118  
CE 1434




**TEST-IT Gel**  
Anti-DVI+  
REF TG119  
CE 1434



**TEST-IT Gel**  
Anti-IgG  
REF TG120  
CE 1434



**TEST-IT Gel**  
Anti-Kell  
REF TG121  
CE 1434



**TEST-IT Gel**  
ABO Type  
REF TG122  
CE

PACKAGING TYPE	SHELF LIFE	STORAGE
40 Cards	15 Months	2 - 25 °C

## REAGENTS AND SOLUTIONS

REFERENCE CODE	PRODUCT NAME	PACKAGING TYPE	SHELF LIFE	STORAGE
TG301	TEST-IT LISS	100 ml	2 Years	2 - 25 °C
TG302	TEST-IT LISS	500 ml	2 Years	2 - 25 °C
TG303	TEST-IT Bromelin	30 ml	9 Months	2 - 8 °C
TG304	TEST-IT Bromelin	100 ml	9 Months	2 - 8 °C
TG305	TEST-IT Papain	10 ml	6 Months	2 - 8 °C



**TEST-IT  
PAPAIN**



**TEST-IT  
BROMELIN**



**TEST-IT  
LISS**

Papain, Bromelin: Liquid Bromelin and Papain solutions for performing enzyme tests using gel techniques.

LISS: A reagent used for preparing red blood cell suspensions used in gel techniques.

## REAGENT RED BLOOD CELLS RBC 0.8%

REFERENCE CODE	PRODUCT NAME	PACKAGING TYPE	SHELF LIFE	STORAGE
TG201	TEST-IT Cell A1-A2-B-O	4 x 10 ml	6 Weeks	2 - 8 °C
TG202	TEST-IT Cell Reverse A1-B	2 x 10 ml	6 Weeks	2 - 8 °C
TG203	TEST-IT Cell Screen 2	2 x 10 ml	6 Weeks	2 - 8 °C
TG204	TEST-IT Cell Screen 3	3 x 10 ml	6 Weeks	2 - 8 °C
TG205	TEST-IT Cell Screen 4	4 x 10 ml	6 Weeks	2 - 8 °C
TG206	TEST-IT Cell Screen 2P	2 x 10 ml	6 Weeks	2 - 8 °C
TG207	TEST-IT Cell Screen 3P	3 x 10 ml	6 Weeks	2 - 8 °C
TG208	TEST-IT Cell Screen 4P	4 x 10 ml	6 Weeks	2 - 8 °C
TG209	TEST-IT Cell Panel 11	11 x 5 ml	6 Weeks	2 - 8 °C
TG210	TEST-IT Cell Panel 11 P	11 x 5 ml	6 Weeks	2 - 8 °C
TG211	TEST-IT Cell Coombs Control	1 x 10 ml	6 Weeks	2 - 8 °C
TG212	TEST-IT Cell Quality Control - 2	2 x 6 ml	6 Weeks	2 - 8 °C
TG213	TEST-IT Cell Quality Control - 4	4 x 6 ml	6 Weeks	2 - 8 °C



**TEST IT CELL  
A1**



**TEST IT CELL  
B**

# MANUAL SYSTEMS

REFERENCE CODE	PRODUCT NAME	PACKAGING TYPE
TG401	TEST-IT Gel Centrifuge	1 x 1
TG402	TEST-IT Gel Incubator	1 x 1
TG403	TEST-IT Gel Reader	1 x 1
TG404	TEST-IT Gel Pipette	1 x 1
TG405	TEST-IT Gel Dispenser	1 x 1
TG406	TEST-IT Gel Workstation	1 x 1
TG407	TEST-IT Gel Pipette Tip	500 pcs/pack
TG408	TEST-IT Gel Test Tube	500 pcs/pack



**TEST-IT GEL WORK STATION**  
Racks for tubes and gel cards.



**TEST-IT GEL CENTRIFUGE**  
Gel card centrifuge for 24 cards.



**TEST-IT GEL INCUBATOR**  
Gel card incubator for 24 cards.



**TEST-IT GEL PIPETTE**  
Repeater Pipette (10µl, 25µl, 50µl)



**DISPENSER**  
500µl

- \* Predefined speed and time : 990 rpm 9 minutes
- \* LCD display shows remaining processing time
- \* Removable and exchangeable spinning head
- \* Electrical lid lock, over speed and imbalance protection.

- \* 2 independent incubation zones (2x12)
- \* Standard program: 15 minutes at 37 °C
- \* Adjustable time
- \* LCD display shows remaining processing time

- Dispensing volumes;
- \* 10 µL to 40 dispensing
  - \* 25 µL to 16 dispensing
  - \* 50 µL to 8 dispensing

- \* Fixed volume 500 µL
- \* Attachable different types of bottles
- \* Easy to use volume control

## TEST-IT GEL READER GEL CARD READER



### Barcode Systems:

- \* Identify TEST-IT Gel Cards by barcode  
Card code, card name and card barcode
- \* External barcode scanner  
Sample barcode and sample name

### Software:

- \* Digitize and process card images
- \* Interpret results automatically

### Data Management:

- \* Record and store data to ensure traceability
- \* Automated upload of results to LIS/HIS.

### User Menu:

- \* Touch Screen
- \* Mini Printer
- \* USB Interface

### HEADQUARTERS

İTOB 10017 Sok. No:2 Menderes-İzmir-TURKIYE

### FACTORY 2:

İTOB 10031 Sok. No:15 Menderes-İzmir-TURKIYE

### FACTORY 3:

İTOB OSB 10001 Sk. No:2 Menderes-İzmir -TURKIYE

### TEST-IT GmbH

Edelweißstraße 4 München/ GERMANY

### PHONE

+90 232 376 80 81

### FAX

+90 232 376 80 40

### E-MAIL

info@turklab.com.tr

**TÜRLAB**

TEST IT



www.turklab.com.tr

