

1 Second Testing

- QUICK
- SIMPLE
- ACCURATE

6 IN 1

- BODY TEMPERATURE
- SURFACE TEMPERATURE
- AMBIENT TEMPERATURE
- HIGH TEMPERATURE ALARM
- MEMORIZE 32 TEST DATA
- THREE COLOR BACKLITS

TECHNICAL SPECIFICATIONS

1. Normal using condition:

Ambient temperature: 10°C ~ 40°C (50°F ~ 104°F)

Relative humidity: ≤ 85%

Pressure altitude: 700hPa ~ 1060hPa

2. Batteries: DC 3V (2 pcs AA batteries)

3. Unit size: 155 x 100 x 40mm (L x W x H)

4. Unit weight (without battery): 105g

5. Temperature display resolution: 0.1°C (0.1°F)

6. Measuring range:

In body mode: 32.0°C ~ 42.9°C (89.6°F ~ 109.2°F)

Under body mode, there is three color backlits:

Green color backlit:

≤ 37.3°C (99.1°F), means normal temperature.

Orange color backlit:

37.4°C ~ 37.9°C (99.3°F ~ 100.2°F), means low fever.

Red color backlit:

≥ 38°C (100.4°F), means high fever.

In surface temp mode: 0°C ~ 60°C (32°F ~ 140°F)

In room mode: 0°C ~ 40°C (32°F ~ 104°F)

7. Precision:

32.0°C ~ 34.9°C (89.6°F ~ 94.8°F) ± 0.3°C (± 0.6°F)

35.0°C ~ 42.0°C (95.0°F ~ 107.6°F) ± 0.2°C (± 0.4°F)

42.1°C ~ 42.9°C (107.8°F ~ 109.2°F) ± 0.3°C (± 0.6°F)

8. Consumption: ≤ 300mW

9. Accuracy: ± 0.2°C (0.4°F)

10. Measuring distance: 3cm ~ 5cm (1.2in ~ 2in)

11. Automatic power-off: < 30 secs

12. Memory: 32 sets

13. Measuring position: Forehead skin



Berrcom®

GUANGZHOU BERRCOM MEDICAL DEVICE CO., LTD.

Address: No.38 Huanzhen Xi Road, Dagang Town, Nansha, 511470 Guangzhou, Guangdong, PEOPLE'S REPUBLIC OF CHINA

Tel: +86(20)34938449

Fax: +86(20)34936960



Please read the product manual carefully or purchase and use under the guidance of medical staff
For details on contraindications or precautions