

Gastroscopy



- Using in diagnostic and surgical procedures in gastroenterology, upper gastrointestinal endoscopy, exploration and treatment of tumors and polyps, etc.

Gastroscopy

Model No.	Distal End Diameter (mm)	Insertion part (mm)	Working channel (mm)	Deflection	Working Length (mm)	Field of View	Depth of Field (mm)
SP-WJA9530	Φ10.5	Φ9.5	3.0	Up 180°Down 90° Right 90°Left 90°	1050	140°	3-100
SP-WJA7020	Φ8.0	Φ7.0	2.0	Four-way 180°			