



*Where*  
Healing Begins

**PhysioPlus**<sup>®</sup>  
REHAB & SPORT

**PRODUCT CATALOGUE**

[www.istanbulsehirmedikal.com](http://www.istanbulsehirmedikal.com)

# Where Healing Begins

PHYSIOPLUS, Mobility Aids, orthopedics designer & manufacturer since 2014.

Our Vision is to be a major player in the world of medical equipment. The PARTNER has specialized over the years in the manufacturing and designing of orthopedics, wheelchairs, Commode chairs, comfort chairs, Crutches, Canes, Rollator and any other product that may accompany a patient or user in maintaining its autonomy and mobility.

Our European development and design teams are equipped with precision high-tech tools resulting in reliability and quality for each product we produce.

Today, PARTNER is a group of more than 100 employees with subsidiaries in Europe and Asia and delivers products sold in more than 15 countries worldwide.

PhysioPlus<sup>®</sup>  
REHAB & SPORT



# C O N T E N T S

**CERVICAL COLLARS**  
01-03

**SHOULDER SUPPORTS**  
04-07

**ELBOW SUPPORTS**  
08-09

**WRIST AND FINGER SPLINTS**  
10-17

**CORSETS**  
18-23

**KNEE SUPPORTS**  
24-31

**ANKLE SUPPORTS**  
32-36

**FOOT CARE**  
37-45

**PEDIATRIC PRODUCTS**  
46-52

**OTHER PRODUCTS**  
53-55





# CERVICAL COLLARS



# CERVICAL COLLARS

## CERVICAL COLLAR (Sponge)

Made of compressed Sponge witch is covered by cotton fabric. Keeps the neck warm and stabilize the neck and the head.

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
After the traumas and bruises.



**PH-1001**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm
XL	40-42 cm



## NELSON COLLAR



Made of the foam witch has holes in it that helps cooling the neck, allows the x-ray to goes throw it, non-allergic, reduces the pain by reducing the ability of the head and the neck to move.

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
After the traumas and bruises.

**PH-1002**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm
XL	40-42 cm



## PHILEDELPHIA COLLAR

Made of the foam, reduces the ability of moving the neck and head thanks to the plastics supports on the front or the back, allows the x-ray to go throw it.

### USAGE

Soft tissue trauma.  
Ligaments and bones injuries.  
Sliding discs.  
Stabilizations in the cases of the uncon-  
sciousness.



**PH-1003**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm



## PHILEDELPHIA COLLAR (Front Open)



Made of the foam, reduces the ability of moving the neck and head thanks to the plastics supports on the front or the back, allows the x-ray to go throw it, the front has a hole in the Trachae area to be used in case of emergency.

### USAGE

Soft tissue trauma.  
Ligaments and bones injuries.  
Sliding discs.  
Stabilizations in the cases of the uncon-  
sciousness.

**PH-1004**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm
XL	40-42 cm



## CERVICAL COLLAR WITH / WITHOUT CHIN SUPPORT (Plastic)

Made of PVC with chin support in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
Cases of radiculopathy  
After the traumas and bruises



**PH-1005 / PH-1006**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm
XL	40-42 cm



## WIRE-FRAME COLLAR JR



Made of the aluminum covered by polyurethane to stabilize the neck as wished

### USAGE

Discs hernia.  
Arthritis.  
Recovery from the surgery.  
Neck injuries.  
Tissue injuries

**PH-1007**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm
XL	40-42 cm



## EMERGENCY NECK COLLAR

Made of the PVC plastic with a foam layer on the inner face of it.  
Can be adjusted to 3 different sizes.  
Has a lock system that helps stabilizing the neck.  
Has a hole in front of it to provide faster access to the neck in cases of emergency.  
Allows the X-ray to go through it.

### USAGE

Stabilizing the neck during the emergencies.  
In coma cases.



**PH-1008**

### Size:

S	34-36 cm
M	36-38 cm



## MINERVA ORTHOSIS



Made of the Thermoplastic with is covered from the inside by sponge that help stabilizing the neck and the upper part of the body.  
Stabilize the chest and reduces the neck's ability of movement.

### USAGE

After the surgeries of the muscles and ligaments.  
Arthritis.  
Bone fractures.

**PH-1009**

### Size:

S	34-36 cm
M	36-38 cm
L	38-40 cm



# SHOULDER SUPPORTS





## ARM SLING (Mesh Fabric Black)

Made of the 100% Polyester, nice on the skin, two layers of fabric provides more stabilization for the shoulder and the elbow and reduces the pressure on them.

**USAGE**  
Shoulder dislocations.  
After surgeries  
prevent edemas  
Fractures



**PH-2001**

Size:

S	34-38 cm
M	38-44 cm
L	44-50 cm
XL	50-58 cm
2XL	58-66 cm
3XL	66-72mm



## ARM SLING LUX



Made of sponge that is covered with the terry fabric from the inside. Stabilize the elbow and the shoulder and reduces the pressure on them.

**USAGE**  
Shoulder dislocations  
After surgeries  
prevent edemas  
Fractures

**PH-2002**

Size:

S	34-38 cm
M	38-44 cm
L	44-50 cm
XL	50-58 cm
2XL	58-66 cm
3XL	66-72mm



## ARM SLING GABARDIAN

Made of Gabardian fabric. Stabilize the elbow and the shoulder and reduces the pressure on them.

**USAGE**  
Fractures and bruises of the arm.  
Bruised shoulder.



**PH-2003**

Size:

S	34-38 cm
M	38-44 cm
L	44-50 cm
XL	50-58 cm
2XL	58-66 cm
3XL	66-72mm



## VELPAU BANDAGE UNIVERSAL



Made of sponge that is covered with the terry fabric from the inside. Stabilize the shoulder and the upper side of the arm. The straps around the chest provides stabilization in case of dislocations

**USAGE**  
Shoulder dislocations  
After arm and shoulder surgeries  
prevent edemas  
Fractures

**PH-2004**

Size:  
**UNISIZE**



## PADDED ARM SLING-(45°-60°)

Made of sponge that is covered with the terry fabric from the inside.  
Stabilize the shoulder and the upper side of the arm on 45°-60D.  
The straps around the chest provides stabilization in case of dislocations

### USAGE

Shoulder dislocations  
After surgeries  
prevent edemas  
Fractures  
Armpits infections



**PH-2005**

**Size:**  
**UNISIZE**

## PADDED ARM SLING-30°



Made of sponge that is covered with the terry fabric from the inside.  
Stabilize the shoulder and the upper side of the arm on 30D.  
The straps around the chest provides stabilization in case of dislocations

### USAGE

After surgeries  
prevent edemas  
Fractures

**PH-2006**

**Size:**  
**UNISIZE**

## CLAVICULA BANDAGE

Made of sponge that is covered with the terry fabric from the inside.  
Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.

### USAGE

Clavicle dislocations  
Clavicle fractures



**PH-2007**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## CLAVICULA BANDAGE



Made of sponge that is covered with the terry fabric from the inside.  
Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.


### USAGE

Clavicle dislocations  
Clavicle fractures

**PH-2008**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## CLAVICULA BANDAGE

Made of sponge that is covered with the terry fabric from the inside. Stabilize the clavicle thanks to the straps that pulls the shoulder to the back to give the appropriate position.

### USAGE

Clavicle dislocations  
Clavicle fractures



### PH-2009

#### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## SHOULDER SUPPORT (Neoprene)



Made of perforated neoprene fabric in 3 mm thickness. Provides warmth to the shoulder with protection in cases of the minor pain.

### USAGE

Shoulder dislocations  
After surgeries  
Warming the Shoulder

### PH-2010

#### Size:

S	34-38 cm
M	38-44 cm
L	44-50 cm
XL	50-58 cm
2XL	58-66 cm



## SHOULDER SUPPORT (Wool)

Made of elastic cotton fabric specially de-signed to reduce sweating. Its inside is covered by smooth non-allergic fabric. Designed to fit under clothes.

### USAGE

Warming the Shoulder.



### PH-2011

#### Size:

S	34-38 cm
M	38-44 cm
L	44-50 cm
XL	50-58 cm
2XL	58-66 cm



## SHOULDER ORTHOSIS



Made of thermoplastic covered by terry fabric from the inside. Has an optional elastic ligament. Stabilize the shoulder and the upper side of the arm.

### USAGE

Stabilizing the shoulder and the upper part of the arm.  
After surgeries.  
Protection and Stabilization.

### PH-2012

#### Size:

S	32-35 cm
M	36-39 cm
L	40-43 cm



# ELBOW SUPPORTS



## TENNIS ELBOW SUPPORT

Made of neoprene fabric that provides warmth and protection. Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epicondyle heads.

### USAGE

Lateral epicondylitis.  
Internal articular fluid. Elbow pain.



**PH-3001**

### Size:

S	22-25 cm
M	25-28 cm
L	28-31 cm
XL	31-34 cm
2XL	34-37 cm



## EPICONDYLITIS BANGAGE UNISIZE

Made of neoprene fabric that provides warmth and protection. Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epicondyle heads.

### USAGE

Used in the treatment of epicondylitis by decreasing the load on the tendons by applying compression.



**PH-3002**

### Size:

**UNISIZE**

## PADDED ELBOW SUPPORT

Made of neoprene fabric that provides warmth and protection. Used by athletes in cases of minor pains.

**USAGE**  
Protects elbow area against impacts.



**PH-3003**

Size:

S	22-25 cm
M	25-28 cm
L	28-31 cm
XL	31-34 cm
2XL	34-37 cm



## TENNIS ELBOW SUPPORT UNISIZE



Made of thermoplastic covered by terry fabric from the inside. The angle can be adjusted easily The Upper part can be adjusted in (0-15-30-45-60-75-90-105) degrees The Lower part can be adjusted in (0-15-30-45-60-75-90-105) degrees

**USAGE**  
Enables controlled mobilization of elbow. medial and lateral epicondyle fractures. After Surgeries. Elbow dislocation, elbow arthrolysis.

**PH-3006**

**Size:**  
**UNISIZE**

## ELBOW CONTRACTURE ORTHOSIS

Made of thermoplastic covered by terry fabric from the inside. The angle can be adjusted easily The Upper part can be adjusted in (0-15-30-45-60-75-90-105) degrees The Lower part can be adjusted in (0-15-30-45-60-75-90-105) degrees

**USAGE**  
Enables controlled mobilization of elbow. medial and lateral epicondyle fractures After Surgeries. Elbow dislocation, elbow arthrolysis.



**PH-3007**

**2 Size (S-M)**  
**2 Side (Right-Left)**

## ELBOW IMMOBILIZER



Made of sponge that is covered with the terry fabric from the inside. Stabilize the Elbow thanks to the metal pieces in it

**USAGE**  
Protections at impacts. Fractures. Protection during training and sport activities.

**PH-3008**

**Size:**  
**UNISIZE**



# WRIST AND FINGER SIPLINTS



## WRIST SPLINT LONG

Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside. Provides better stabilization thanks to the elastic bandage on the wrist.

**USAGE**  
Provides immobilization and constrains to the movements of the wrist  
Arthritis.  
Rheumatism.



### PH-4001

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## WRIST SPLINT (NEOPRENE) LONG



Made of neoprene fabric that provides warmth and protection. Constrains movements of the wrist thanks to the aluminum support inside.

**USAGE**  
Wrist tunnel cases.  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.

### PH-4002

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## WRIST SPLINT (NEOPRENE) LONG Right and Left

Made of neoprene fabric that provides warmth and protection. Constrains movements of the wrist thanks to the aluminum support inside.

**USAGE**  
Wrist tunnel cases.  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.



### PH-4004

Size:

**UNISIZE**

## WRIST SPLINT WITH THUMB SUPPORT LONG



Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside. The thumb supporting rod stabilizes it in the rest position. Provides better stabilization thanks to the elastic bandage on the wrist.

**USAGE**  
Protection and stabilization of the wrist and the thumb.  
Arthritis and rheumatism.  
Wrist tunnel cases.  
Pains in the thumb and wrist.

### PH-4005

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## WRIST SPLINT (NEOPRENE) With Thumb Support Long

Made of neoprene fabric that provides warmth and protection. It has an aluminum support with a bandage to spread the pressure all over the hand.

### USAGE

Stabilizing the thumb.  
Wrist tunnel cases.  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.



**PH-4006**

### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## WRIST SPLINT SHORT



Made of elastic cotton fabric specially designed to constrain the movements of the wrist thanks to the aluminum rod inside.

### USAGE

Protection and stabilization of the wrist  
Arthritis and rheumatism.  
Wrist tunnel cases.  
Pains in wrist.

**PH-4007**

### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## WRIST SPLINT SHORT NEOPRENE

Made of neoprene fabric that provides warmth and protection. It has an aluminum support to spread the pressure all over the hand.

### USAGE

Protection and stabilization of the wrist  
Arthritis and rheumatism.  
Wrist tunnel cases.  
Pains in wrist



**PH-4008**

### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## WRIST SPLINT WITH THUMB SUPPORT (Short) Neoprene Fabric



Stabilizing the thumb's joint. Thumb ligament's problems.

### USAGE

Stabilizing the thumb's joint.  
Thumb ligament's problems.

**PH-4009**

### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm
XL	21-23 cm
2XL	23-25 cm



## ADJUSTABLE WRIST SPLINT (NEOPRENE)

It is manufactured from perforated fabric. It is a type of splint that provides immobilization of the wrist in neutral position. It has an aluminium support inside and 2 plastic supports at the top. There is a bandage that provides stabilization in the laced system around the wrist and can be used in both hands.

### USAGE

Carpal Tunnel Syndrome  
In tendon and soft tissue injuries  
Rehabilitation after surgical intervention



**PH-4011**

**Size:**  
**UNISIZE**

## NEOPRENE WRIST BANDAGE



It is a simple wrist bandage. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

### USAGE

In wrist strains.  
To prevent sports injuries.  
In carpal dislocation.

**PH-4013**

**Size:**  
**UNISIZE**

## NEOPRENE HAND WRIST SUPPORT

It is a simple wrist bandage that put thumb to inside. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

### USAGE

Twists, Trumas, bruises,  
Protection while doing sports.  
Protection at work when the jobs involves using the wrist a lot.



**PH-4014**

**Size:**  
**UNISIZE**

## THUMB SUPPORT (DE QUEUVAIN)



It is a type of splint that immobilize thumb. The plastic stele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

### USAGE

CMC and MCP joint irritation.  
In de quervian syndrome.  
In ligament damage.

**PH-4015**

**Size:**  
**UNISIZE**



## THUMB SUPPORT

Made of thermoplastic material, it is a type of splint that immobilize thumb. The plastic. Male inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

### USAGE

CMC and MCP joint irritation.  
In de quervain syndrome.  
In ligament damage.



**PH-4016**

### Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



## WRIST BAND WITH THUMB SUPPORT



It is a type of splint that immobilize thumb, The aluminium atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

### USAGE

CMC and MOP joint irritation.  
In de quervain syndrome.  
In ligament damage.

**PH-4017**

**Size:**  
**UNISIZE**

## WRIST BAND WITH THUMB SUPPORT LONG

It is a type of splint that immobilize thumb. The aluminium atele inside provides stabilization of the thumb. It can adjust with velcro on it for wrist.

### USAGE

CMC and MCP joint irritation.  
In de quervain syndrome.  
In ligament damage.



**PH-4018**

**Size:**  
**UNISIZE**

## WRIST SPLINT WITH DOUBLE SIDE



There is an anatomically shaped aluminium support that provides immobilization of the wrist in the neutral position. There is a bandage that increases the stabilization around the wrist. It is a standard size and can be personally adjusted with velcro.

### USAGE

Carpal Tunnel Syndrome.  
In tendon and soft tissue injuries.  
Rehabilitation after surgical intervention.

**PH-4019**

**Size:**  
**UNISIZE**

## WRIST SPLINT WITH THUMB SUPPORT NEOPRENE

There is an aluminium support under the wrist and near the thumb. It can use to both hands by changing aluminium support under the wrist. There is a velcro that increases the stabilization around the wrist.

### USAGE

In de guervain syndrome  
Carpal Tunnel Syndrome  
Rehabilitation after surgical intervention  
In wrist and thumb calcification  
In tendon and soft tissue injuries



**PH-4020**

**Size:**  
**UNISIZE**

## THERMOPLASTIC RELIEVING SPLINT



Wand and wrist orthosis manufactured from thermoplastic material. Keeps the fingers and the wrist in recovery position. Enables immobilization of wrist, used for resting position following plaster treatment of fractures and painful positions of the wrist.

### USAGE

For bone fractures of the wrist after plaster casting  
For resting in soft tissue injuries  
Carpal Tunnel Syndrome  
Rehabilitation after surgical intervention  
In tendon and soft tissue injuries

**PH-4021**

**Size:**  
S 15-17 cm  
M 17-19 cm  
L 19-21 cm



## THERMOPLASTIC ANTI-SPASTICITY HAND SPLINT

Made of thermoplastic covered by terry fabric from the inside. Stabilizes the wrist and the fingers in a relaxing opened position.

### USAGE

Spasticity  
Wrist problems,



**PH-4022**

**Size:**  
S 15-17 cm  
M 17-19 cm  
L 19-21 cm



## THERMOPLASTIC ANTI-SPASTICITY HAND SPLINT LUX



Made of thermoplastic covered fabric from the inside. Stabilizes the wrist and the fingers in a relaxing opened position.

### USAGE

Spasticity  
Wrist problems.

**PH-4023**

**Size:**  
S 15-17 cm  
M 17-19 cm  
L 19-21 cm



## THERMOPLASTIC DYNAMIC HAND-FINGER SPLINT

Made of thermoplastic covered by terry fabric from the inside, Pulls up the fingers and let them move a bit to help training them.

**USAGE**  
Rehabilitation.  
Extensor tendon.s.



**PH-4024**

Size:

S	15-17 cm
M	17-19 cm
L	19-21 cm



## THERMOPLASTIC KLEINERT HAND-FINGER SPLINT



Made of thermoplastic covered by terry fabric from the inside. Pulls up the fingers and let them move a bit to help training them.

**USAGE**  
Rehabilitation.  
Flexure tendon.

**PH-4025**

Size:

S	15-17 cm
M	17-19 cm



## MALLET FINGER SPLINT

Made of high quality plastic. Stabilize the upper part of the finger and prevent it from moving.

**USAGE**  
Mallet finger deformities.  
Interruption of tendons.  
Fractures and bruises.



**PH-5001**

Size:

S	4.5 - 5 cm
M	5 - 6 cm
L	6 - 7 cm
XL	7 - 8 cm
2XL	8 - 9 cm

## FROG SPLINT



Made of the Soft aluminum witch has the ability to reshape as wished. The inner side is covered with the plastoz witch helps stabilizing the finger.

**USAGE**  
Finger fractures.  
Strings and nerve injuries.  
Finger injuries and pain.

**PH-5002**

Size:

S	4.5 - 5 cm
M	5 - 6 cm
L	6 - 7 cm
XL	7 - 8 cm
2XL	8 - 9 cm

## BASEBALL SPLINT

Made of the Soft aluminum witch has the ability 10 reshape as wished. The inner side is covered with the plastoz witch helps stabilizing the finger.

### USAGE

Finger fractures,  
Strings and nerve injuries.  
Finger injuries and pain.



**PH-5003**

### Size:

S 4 - 6 cm

M 6 - 8 cm

L 8 - 10 cm

## BUTONIER FINGER SPLINT

Made of The Soft aluminium witch has the ability to reshap as wished. The inner side is covered with the plastoz witch helps stabilizing the finger.

### USAGE

Rehabilitation.  
Strings injuries.  
Finger movement training.



**PH-5004**

### Size:

S 4 - 6 cm

M 6 - 8 cm

L 8 - 10 cm

## STRAIGHT FINGER SPLINT

Made of the Soft aluminium which has the ability to reshape as wished, The inner side is covered with the plastoz witch helps stabilizing the finger.

### USAGE

Finger fractures.  
Strings and nerve injuries.  
Finger injuries and pain.

RESMİ  
YOK

**PH-5005**

### Size:

Standard



# CORSETS



## ELASTIC LUMBOSTAD CORSE-WITH BELT (26 CM)

Made of elastic cotton fabric specially designed to reduce sweating. Has 2 metal and 2 plastic supportors reduces the pressure on the discs. Has elastic belt provides extra stabilization, Its 26 cm High. The metal supportors can be removed while washing.

### USAGE

Lumbar problems.  
Chronic back pain.



**PH-6001**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## LUMBOSACRAL CORSET 26 CM



Made of elastic cotton fabric specially designed to reduce sweating. Has 2 metal and 2 plastic supportors reduces the pressure on the discs. It's 26 cm High. The metal supportors can be removed while washing.

### USAGE

Sliding discs.  
Chronic back pain.  
Before and after surgeries,  
Lumbar spine pain.

**PH-6002**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## ELASTIC LUMBOSTAD CORSE-WITH BELT (32 CM)

Made of elastic cotton fabric specially designed to reduce sweating. Has 2 metals and 2 plastics supportors reduces the pressure on the discs. Has elastic belt provides extra stabilization. It's 32 cm High. The metal supportors can be removed while washing.

### USAGE

Sliding lumbar and sacral discs.  
Chronic back pain. Lumbar problems. Before and after surgeries.  
Reduces the pressure on the discs.



**PH-6003**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## DORSO LUMBAR CORSET THORACO LUMBAR)



Made of elastic cotton fabric specially designed to reduce sweating. Has 2 metal and 2 plastic supportors at the lower back. Supports the spine and the shoulders thanks to its design. Comfortable during daily use.

### USAGE

Spine problems.  
Osteoporosis.  
The small humping.  
Support before and after Surgeries.

**PH-6004**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



# CORSETS

## POSTUREX BANDAGE

Made of elastic cotton fabric specially designed to reduce the Sweating.  
Has a plastic supporting reduces the ability to bend.  
Has side laces helps pulling the shoulders to keep the body's stabilization.  
it can be easily worn under the clothes.

### USAGE

Shoulder and back muscles pain.  
Other muscle problems,  
Long computer and car driving hours.  
Back bags carrying.



**PH-6005**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## MAGNETIC BACK CORSET



Adjustable straps design to fit all sizes the flexibility that gives you the freedom of movement. It has a magnetic belt that supports the spine to provide stability for those suffering from lower back pain.

### USAGE

Great for lower back pain and supporting the back. It also keeps lower back in an upright or extended position it reduces pain in muscles, ligaments, discs and nerve roots of the lower back.

**PH-6006**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## ABDOMINAL CORSET

Made of elastic cotton fabric specially designed to reduce sweating  
Stabilize the abdomen.  
Can be adjusted as needed,  
Reduces the the pains in the abdomen area. Designed to fit under clothes.

### USAGE

Post abdominal surgery and birth  
Peritoneal injury  
Abdominal sagging  
Post umbilical and ventral hernia surgery




**PH-6007**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## CHEST CORSET



Made of elastic cotton fabric specially designed to reduce sweating.  
Its inside is covered by smooth nonallergenic fabric.  
Designed to fit under clothes.

### USAGE

Sternal Support after coronary surgery  
Costa and sternum fractures  
Post chest surgery and trauma

**PH-6008**

### Size:

S	80-90 cm
M	90-100 cm
L	100-110 cm
XL	110-120 cm
2XL	120-130 cm





## NEOPRENE LUMBOSACRAL CORSET WITH BELT

Made of neoprene fabric that provides warmth and protection.  
Has 2 metals and 2 plastics supporters reduces the pressure on the discs.  
Has a belt provide extra warmth and pressure.

### USAGE

Sliding lumbar and sacral discs,  
Back pain.  
Prevent edemas.  
Enhance the blood circulation.



**PH-6009**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## NEOPRENE LUMBOSACRAL CORSET 26 CM



Made of neoprene fabric that provides warmth and protection. Has 2 metals and 2 plastics supporters reduces the pressure on the discs. Its 26 cm High.  
Has a belt provide extra warmth and pressure.

### USAGE

Sliding discs.  
Chronic back pain.  
Before and after surgeries.  
Lumbar spine pain.

**PH-6010**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## ABDOMINAL CORSET (NEOPRENE)

Made of neoprene fabric that provides warmth and protection.

### USAGE

Provides warmth to the abdomen and the lower part of the back.  
Helps burning the fats and losing weight during doing sports.



**PH-6011**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## HYPEREXTENSION CORSET



Made of a special designed aluminum. The inner side is covered with a smooth non-allergic sponge fabric. The 3 points system stabilize the spine and prevent its movement which speeds up the healing process.

### USAGE

Stable compression fracture of lumbar and lower Thoracic vertebrae  
After inter-vertebral disc surgery  
After non-stable vertebral body fracture surgery. Non-operable spinal tumor and its metastases. Post-operation orthopedic care and orthosis  
After spinal decompression surgery

**PH-6012**

### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm





# CORSETS

## HYPEREXTENSION CORSET ADJUSTABLE

Made of a special designed aluminum, The inner side is covered with a smooth non-allergic sponge fabric. The special design helps modifying it as needed. Stabilizes the spine and keeps it straight.

### USAGE

Stable compression fracture of lumbar and lower thoracic vertebrae  
After inter-vertebral disc surgery  
After non-stable vertebral body fracture surgery. Non-operable spinal tumor and its metastases. Post-operation orthopedic care and orthosis  
After spinal decompression surgery.



### PH-6013

#### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## COLOSTOMY CORSET



Made of elastic cotton fabric specially designed to reduce sweating. Has a side hole in it so the colostomy base can be attached. Designed to fit under clothes.

### USAGE

Protecting the ostomy collection bag  
Post-ostomy rehabilitation and support purposes  
Decrease risks of parastomal hernia.

### PH-6014

#### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## ABDOMINAL CORSET (WOOL)

Made of 2 layer wool fabric.

### USAGE

Provides warmth to the abdomen and the lower part of the back. Cal and ventral hernia surgery.



### PH-6015

#### Size:

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



## WOOL CORSET (OPEN FRONT SIDE)



Made of 2 layer wool fabric.


### USAGE

Provides warmth to the abdomen and the lower part of the back.

### PH-6016

#### Size:

S	80-90 cm
M	90-100 cm
L	100-110 cm
XL	110-120 cm
2XL	120-130 cm



## INGUINAL RAPTURE BELT RIGHT OR LEFT

Made of a special cotton fabric. Does not effect the daily life. Provides support after surgeries. Fits perfectly under clothes.

**USAGE**

Single sided inguinal hernia.



**PH-6018**

**Size:**  
**Standard**

## INGUINAL RAPTURE BELT (DOUBLE)

Made of a special cotton fabric. Does not effect the daily life. Provides support after surgeries. Fitts perfectly under clothes.

**USAGE**

Double sided inguinal hernia.



**PH-6019**

**Size:**  
**Standard**

## (NEOPRENE) SLIMMING SHORT

Made of neoprene fabric contains pores that help for ventilation effectively improve blood circulation. Keeps the heat, reduces the pressure, reduces the severity of muscle's pain, Helps to slim the hip, thigh and torso. Designed for comfort, support and protection for hamstring and thigh injury. Zipper provide ease at wear.

**USAGE**

Providing heat and compression Inguinal problems. Slimming around the waist and hip line injuries to the thigh.



**PH-6020**

**Size:**

S	70-80 cm
M	80-90 cm
L	90-100 cm
XL	100-110 cm
2XL	110-120 cm



# KNEE SUPPORTS



## (NEOPRENE) KNEE SUPPORT CLOSED PATELLA

Made of neoprene fabric  
Contains pores that help for ventilation, Support knee joint, reduce pain by relieving the pressure on the joint of the knee. Sustains body's own release of heat effectively improve blood circulation, promote metabolism, enhance immunity also applies compression, prevents edema.

**USAGE**  
Arthritis. Intra articular liquid increase Pain induced by knee cap Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities When heat and compression is desired



**PH-7001**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

## (NEOPRENE) KNEE SUPPORT CLOSED WITH SILICONE PAD



Made of neoprene fabric  
Contains pores that help for ventilation, Support knee joint, reduce pain by relieving the pressure on the joint of the knee. The perimeter of the patella is closed and supported with a silicone pad. Sustains body's own release of heat effectively improve blood circulation, promote metabolism, enhance immunity also applies compression, prevents edema.

**USAGE**  
Arthritis. Intra articular liquid increase Pain induced by knee cap Post-op or trauma for protection and rehabilitation purposes, Protection during and after the sportive activities When heat and compression is desired

**PH-7002**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

## (NEOPRENE) KNEE SUPPORT CLOSED PATELLA

Made of neoprene fabric  
Contains pores that help for ventilation, Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap. Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

**USAGE**  
Arthritis. Intra articular liquid increase Pain induced by knee cap Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities When heat and compression is desired



**PH-7004**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm

## (NEOPRENE) KNEE SUPPORT PATELLA WITH LIGAMENT SUPPORT-SHORT



Made of neoprene fabric  
Contains pores that help for ventilation. Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

**USAGE**  
Arthritis. bursitis, tendinitis  
Sprains and stains  
Meniscus problems  
Post-op or trauma for protection purposes  
Ligament instability  
Sporting injuries

**PH-7004**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## (NEOPRENE) KNEE SUPPORT PATELLA WITH LIGAMENT SUPPORT-LONG

Made of neoprene fabric  
Contains pores that help for ventilation,  
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee,

**USAGE**  
Arthritis, bursitis, tendinitis  
Sprains and stains. Meniscus problems  
Post-op or trauma for protection purposes  
Ligament instability  
Sporting injuries



**PH-7005**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## (NEOPRENE) KNEE SUPPORT HALF OPENED PATELLA WITH LIGAMENT SUPPORT



Made of neoprene fabric  
Contains pores that help for ventilation.  
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Flexible wire on the sides increase stabilization and decrease the pressure on the knee.

**USAGE**  
Arthritis, bursitis, tendinitis  
Sprains and stains. Meniscus problems  
Post-op or trauma for protection purposes  
Ligament instability. Sporting injuries

**PH-7006**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## KNEE SUPPORT WITH HINGED STABILIZING JOINTS

Contains pores that help for ventilation.  
Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension. Velcro fasteners on the upper and lower calves increase stabilization.

**USAGE**  
Ligament injury and instability  
Meniscus lesions  
1 and 2 degree sprains and stains  
Post-op or trauma for protection purposes



**PH-7007**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## HALF OPENED KNEE SUPPORT WITH HINGED STABILIZING JOINTS



Made of neoprene fabric  
Contains pores that help for ventilation. Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension. Velcro fasteners on the upper and lower calves increase stabilization.

**USAGE**  
Ligament injury and instability  
Meniscus lesions  
1 and 2 degree sprains and stains  
Post-op or trauma for protection purposes

**PH-7008**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## KNEE SUPPORT WITH HINGED STABILIZING JOINTS-FULL OPEN FRONT

Made of neoprene fabric  
Contains pores that help for ventilation.  
Aluminum joints on both sides provide maximum support, allow controlling the motion and prevent knee hyperextension.  
Velcro fasteners on the upper and lower calves increase stabilization.  
Velcro fasteners on the upper and lower calves increase stabilization.  
Since it can be completely open from the front side, provides ease of use to the elderly, heavy weight and post-op patients.

**USAGE**  
Ligament injury and instability  
Meniscus lesions  
ve 2 degree sprains and stains  
Post-op or trauma for protection purposes



**PH-7009**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## KNEE SUPPORT PATELLA WITH CROSSED LIGAMENT SUPPORT-NEOPRENE



Made of neoprene fabric  
Contains pores that help for ventilation,  
Contains opened the patella supported by a pad. Medial (inside) and lateral (outside)  
Flexible wire offer overall knee joint stability.  
Transverse straps adhering with 3-point principle provide stabilization far ant. Cruciate ligaments and patella.

**USAGE**  
Ant. Cruciate ligament injury and knee instability  
Patellar stabilization  
Meniscus lesions  
Post-op or trauma for protection and rehabilitation purposes

**PH-7010**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## KNEE SUPPORT FOR PATELLAR TENDON

Made of neoprene fabric that provides warmth and protection.  
The special pad supports the patellar tendon and with the applied compression reduces pain.

**USAGE**  
Os-good Schlatter's disease Patellar tendinitis  
Patellar chondromalacia Patellofemoral pain



**PH-7011**

Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## PATELLAR TENDON BAND UNISIZE



Made of neoprene fabric that provides warmth and protection.  
The special pad supports the patellar tendon and with the applied compression reduces pain.  
Does not restrict joint motion.  
Providing comfortable all-day wear.  
The knee support strap is made with neoprene so it provides a comfortable and flexible fit.  
Used for both knees.

**USAGE**  
Os-good Schlatter's disease  
Patellar tendinitis  
Patellar chondromalacia  
Patellofemoral pain

**PH-7012**

Size:  
Standard

## WOOL KNEE SUPPORT

Manufactured of wool for warmth and compression

### USAGE

Knee pain  
Liming  
Rheumatic pains  
Temperature for



**PH-7013**

### Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## KNEE SUPPORT OPEN PATELLA



Made of neoprene fabric  
Contains pores that help for ventilation.  
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap.  
Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

### USAGE

Arthritis. Intra articular liquid increase Pain induced by knee cap  
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities. When heat and compression is desiredg injuries

**PH-7014**

**Size:**  
**Standard**

## KNEE SUPPORT OPEN PATELLA WITH LIGAMENT

Made of neoprene fabric  
Contains pores that help for ventilation.  
Contains opened the patella supported by a pad. Patellar Knee Support is ideal for patients who has pain on patella and has post-operative pain. Can be used daily as it does not compress your knee cap.  
Recommended for people who have pain in their knee cap with and difficulty when climbing stairs or squatting.

### USAGE

Arthritis. Intra articular liquid increase Pain induced by knee cap  
Post-op or trauma for protection and rehabilitation purposes. Protection during and after the sportive activities. When heat and compression is desiredg injuries



**PH-7015**

**Size:**  
**Standard**

## ADJUSTABLE ANGLE KNEE ORTHESIS WIT BAR LUX



Use for protected activity during post-op rehab following ligamentous surgical repairs and reconstruction. Controlled R.O.M. for post-traumatic, post-surgical repair of human knee. It can be extended and shortened with the buttons on it.

### USAGE

Controlled and gradual joint mobilization after ligament injury operations  
Meniscal repairs. For knee rest and stabilize at the desired angle after surgery.  
In knee recurvatum injuries  
Stable femoral fractures

**PH-7020**

**Size:**  
**Standard**

## ADJUSTABLE ANGLE KNEE ORTHESIS WITH BAR

Adjustable-degree knee orthosis, manufactured from aluminium bar and sponge laminated cloth with a terry inner side. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°, 0°, 15°, 30°, 45°, 60°, 75°, 90°) angle adjustment possible. Provides exceptional support and stability following knee injury. Use for protected activity during post-op rehab following ligamentous surgical repairs and re-construction. Controlled R.D.M. for post-traumatic, post-surgical repair of human knee.

### USAGE

Controlled and gradual joint mobilization after ligament injury operations. Meniscal repairs  
For knee rest and stabilize at the desired angle after surgery.  
In knee recurvatum injuries  
Stable femoral fractures



**PH-7021**

**Size:**  
**Standard (47/58 cm)**

## FUNCTIONAL KNEE BRACE

It is a functional rigid brace that is made of light-weight aluminium material designed according to the anatomical knee shape, providing flexion and extension limitation, Flexion (0°, 15°, 30°, 60°, 90°, 120°) and extension (0°, 15°, 30°) degrees can be adjusted. It provides the most favorable stabilization with the principle of 4 points, It prevents forward and back tibial displacement and prevents knee hyperextension.

### USAGE

Instability caused by collateral ligament rupture  
Conservative treatment of combined instabilities  
Isolated rupture of ant. cruciate ligament or combined ant. and post. Cruciate ligament ruptures



**PH-7022**

**Size:**  
**Standard**

## NEOPRENE KNEE SUPPORT WITH CONTROLLED MECHANISM

Neoprene fabric which Provides firm support and promotes flexibility and it also offers therapeutic heat to soothe aching joints. Controls extension and flexion degrees of the knee movement. Screws that are adjusted in Flexion in degree. 0°, 15°, 30°, 60°, 90°, 120°) and Extension (0°, 5°, 15°, 30°) provides external support by protecting knee junction. Better control, comfort, and increased compliance. And reduces the pain by relieving the load on the joint.

### USAGE

Ligament injuries and instability, Meniscus lesions, 1st and 2nd degree sprains and strains.  
Post-op or trauma for protection and rehabilitation purposes. Degenerative joint disease where extension/flexion control is required.



**PH-7023**

### Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm





## KNEE IMMOBILIZER

Manufactured of sponge laminated with a terry cloth inner coating, This orthosis immobilizes the leg at 20°. Aluminum supports provide immobilization in extension, stabilization in the posterior, medial and lateral. Knee immobilizers are ideal for knee injuries and as a postoperative aid as they are easily removed. Standard size.

Lengths available at 45cm - 55cm and 65 cm.

Thanks to the three-piece pad construction, this orthosis can be easily adjusted to any leg circumference.

### USAGE

Recovery after ligament injury

Meniscus Post-op recovery

Knee joint in trauma cases

Osteochondral fractures

Acute sprains/strains of the knee

Patellar dislocations

Post-op recovery, protection and immobilization



**PH-7024**

**Size:**  
**45 cm**



**PH-7025**

**Size:**  
**45 cm**



**PH-7026**

**Size:**  
**45 cm**

## (NEOPRENE) KNEE SUPPORT CLOSED PATELLA

Made of neoprene fabric  
Contains pores that help for ventilation,  
It supports thigh's muscles  
effectively improve blood circulation,  
Keeps the heat, reduces the pressure.  
reduces the severity of muscle's pain

### USAGE

Muscle injuries to the thigh  
Thigh pain  
Protection during and after the sportive  
activities  
When heat and compression is needed



**PH-7027**

### Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## (NEOPRENE) KNEE SUPPORT CLOSED WITH SILICONE PAD



Made of neoprene fabric  
Contains pores that help for ventilation. It  
supports thigh's muscles  
effectively improve blood circulation, Keeps  
the heat, reduces the pressure, reduces the  
severity of muscle's pain

### USAGE

Muscle injuries to the calf  
Calf pain  
Protection during and after the sportive  
activities  
When heat and compression is needed

**PH-7028**

### Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## (NEOPRENE) KNEE SUPPORT CLOSED PATELLA

It is manufactured to support thigh. It  
reduces pain by radiate heat to thigh area  
Velcros on it provides easy usage. Include  
belt and so it provides comfortable usage to  
patients for daily usage.

### USAGE

Muscle injuries to the thigh  
Thigh pain  
Protection during and after the sportive  
activities  
When heat and compression is needed



**PH-7029**

### Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



## (NEOPRENE) KNEE SUPPORT PATELLA WITH LIGAMENT SUPPORT-SHORT



Hip abduction Orthosis anatomically de-  
signed and manufactured of thermoplastic  
material, providing hip stabilization.  
Parts in contact with the skin are coated I  
with terry cloth. Flexion adjustment possible  
at every 15 degrees, starting from 15 up to 120  
degrees.  
Abduction angle is adjustable.

### USAGE

Stabilization of total hip prosthesis or revision  
Femoral head resection  
Femoral or acetabular osteotomy Recurrent  
hip dislocation

**PH-7030**

### Size:

S	34-36 cm
M	36-38 cm
L	38-41 cm
XL	41-44 cm
XXL	44-47 cm



# ANKLE SUPPORTS



## ANKLE SUPPORT (CROSS STRAP)

Made of neoprene fabric  
Contains pores that help for ventilation,  
The figure 8 elastic transverse strap  
provides extra stabilization and  
compression. The soft neoprene blend  
retains body heat to improve circulation  
and promote healing process.

### USAGE

1st degree twist and torus  
Minor fractures  
Compression on edema  
Protection during and after the sportive  
activities



**PH-8001**

### Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	37-37 cm



## ANKLE SUPPORT WITH SILICONE PAD (MALLEOLUS)



Made of neoprene fabric  
Contains pores that help for ventilation, The  
Crossed lacing provides extra pressure and  
stabilization thanks to the rubber pads on  
the sides.

### USAGE

Ligament instability and rapture Edema  
prevention  
Post-op or trauma for protection and  
rehabilitation purposes

**PH-8002**

### Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	37-37 cm



## LIGAMENT ANKLE SUPPORT

Made of neoprene fabric  
Contains pores that help for ventilation. The  
aluminium around the ankle provides high  
leveled stabilization.  
The aluminium can be modified as wanted.  
The Crossed lacing provides extra pressure  
and stabilization.

### USAGE

Ligament instability and rapture Edema  
prevention  
Post-op or trauma for protection and  
rehabilitation purposes



**PH-8003**

### Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	37-37 cm



## FIGURE OF EIGHT BANDAGE



Made of neoprene fabric  
Contains pores that help for ventilation.  
Contains opened the patella supported by a  
pad. Patellar Knee Support is ideal for  
patients who has pain on patella and has  
post-operative pain. Flexible wire on the  
sides increase stabilization and decrease  
the pressure on the knee.

### USAGE

Arthritis, bursitis, tendinitis  
Sprains and stains  
Meniscus problems  
Post-op or trauma for protection purposes  
Ligament instability  
Sporting injuries

**PH-8004**

**Size:  
Standard**



## ANKLE SUPPORT (PLASTIC + LACED)

Made of neoprene fabric  
Contains pores that help for ventilation.  
Thanks to two specially designed plastic plates located on the malleolus, it provides strong support to the ankle.

### USAGE

Provides stronger stabilization to the ankle with serious ankle sprains.  
Also effective with habitual ankle sprains, chronic ankle pains and can be used as a protector with contact sports.



**PH-8005**

### Size:

S	24-26 cm
M	26-28 cm
L	28-31 cm
XL	31-34 cm
XXL	37-37 cm



## ANKLE SUPPORT (CROSS STRAP)



Made of neoprene fabric  
Contains pores that help for ventilation. The figure 8 elastic transverse strap provides extra stabilization and compression. The soft neoprene blend retains body heat to improve circulation and promote healing process.

### USAGE

1st degree twist and twos  
Minor fractures  
Compression on edema  
Protection during and after the sportive activities.

**PH-8006**

**Size:**  
**Standard**

## ANKLE BRACE WITH AIR PAD

Made with specially designed plastic molds.  
Thanks to the air pads inside, it fits to the foot anatomically and provides maximum stabilization.

### USAGE

Functions as a stabilizer in cases of serious ankle sprains and chronic ankle pains and as a protector in contact sports.  
Also applicable in cases of ankle fractures and in the period after plastering process.



**PH-8007**

**Size:**  
**Standard**

## ANKLE BRACE WITH PLASTIC PAD



The hard outer shell has been made with specially designed plastic molds.  
Thanks to the soft pads inside, it fits to the foot anatomically and provides maximum stabilization.

### USAGE

With its rigid exterior, provides good stabilization in pronation and supination actions. Provides full support for ankle ligament injuries.

**PH-8008**

**Size:**  
**Standard**

## ANKLE SUPPORT (CROSS STRAP)

It comes with 3 pads (12.5°) so it can be adjusted as needed.  
Physiological walk due to special sole.  
Plantar flexion: 0°, 10°, 20°, 30°, 40°  
Dorsiflexion: 20°, 10°, 0°

### USAGE

After ligament, soft tissue and tendon surgeries.  
For post operative and follow up treatments of the talus, calcaneus and medical malleolus fractures.  
For early functional therapy after conservative or surgical treatment for rupture for the Achilles tendon



**PH-8009**

### Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



## ANKLE SUPPORT WITH SILICONE PAD (MALLEOLUS)



It comes with 3 pads (12.5°) so it can be adjusted as needed,  
Physiological walk due to special sole.

### USAGE

After ligament, soft tissue and tendon surgeries.  
For post operative and follow up treatments of the talus, calcaneus and medical malleolus fractures.  
For early functional therapy after conservative or surgical treatment for rupture for the Achilles tendon

**PH-8010**

### Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



## LIGAMENT ANKLE SUPPORT

Ankle orthosis made of plastic material.  
Provides static dorsiflexion assistance and lateral stability for entire anklegoot area.  
Provides ease of use with velcro fasteners going around the ankle.  
Use a scissor to trim insole for exact shoe fit.  
May be modified with a heat gun.

### USAGE

Foot-drop deformity caused by nerve injuries such as hemiplegia, lornber disc hernia  
When 90 degrees ankle angle is needed



**PH-8011**

### Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



## FIGURE OF EIGHT BANDAGE



T shaped ankle orthosis made of thermo-plastic material with terry cloth coating.  
Ensures the fix 90 degrees stand, prevent rotation of the ankle with T Shape,

### USAGE

Stabilizing the position of the leg on hospital patients

**PH-8012**

### Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



## HALLUKS VALGUS NIGHT SPLINT

Foam laminated fabrics, the inner surface is made of terry cloth. Position of the thumb to the desired position is hidden inside the fabric with aluminum splint can be adjusted. designed to be used at Night (in sleep).

**USAGE**  
Helps to keep the big toes aligned, designed to use during night time



**PH-8013**

**Size:**  
**Standard**

## HALLUKS VALGUS DAY SPLINT

Velcro and elastic bandage is made of neo-prene fabric. Day (shoes, slippers ... in) designed to be used at night

**USAGE**  
Helps to keep the big toes aligned, designed to use during day time.



**PH-8014**

**Size:**  
**Standard**

## CASTHOE

with adjustable Velcro tapes can also be applied to both feet.

**USAGE**  
Used instead of shoes in case of plaster feet



**PH-8015**

**Size:**  
S 34-37 cm  
M 38-41 cm  
L 42-45 cm



# FOOT CARE



## SILICONE INSOLE

Made of high quality of silicone material, Specially designed to fit figure of the foot. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Distributes the standing body weight evenly and reduces load and stress on the joints. Absorbs the reaction force of the floor.

### USAGE

Pes planus  
Plantar fasciitis  
Heel pain  
Metatarsal pain



**PH-8017**

### Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



## SILICONE INSOLE WITH ARCH SUPPORT

Specially designed to fit figure of the foot. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.

### USAGE

Provides anatomical support to the sole and enables even distribution of pressure on the sole. Has a relaxing and comforting effect.



**PH-8021**

### Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm





## PERFORATED SILICONE INSOLE WITH ARCH SUPPORT

Specially designed to fit figure of the foot. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.

### USAGE

Provides anatomical support to the sole and enables even distribution of pressure on the sole. Has a relaxing and comforting effect.



**PH-8022**

### Size:

S	34-37 cm
M	38-41 cm
L	42-45 cm



## COVERED SILICONE INSOLES



Specially designed to fit figure of the foot and covered with fabric. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Distributes the standing body weight evenly and reduces load and stress on the joints. Absorbs the reaction force of the floor.

### USAGE

Pee planus. Plantar fasciitis Heel pain. Metatarsal pain

**PH-8023**

### Size:

S	34-37 cm
M	38-39 cm
L	40-41 cm



## CORK INSOLES

Longitudinal and metatarsal arch support manufactured from mix-ture of special cork and rubber material, covered with leather surface. Distributes the standing body weight evenly and reduces load and stress on the joints.



**PH-8032**

### Size:

From	35 to 44
------	----------



## ARC REINFORCED RUBBER INSOLE

Specially designed to fit figure of the foot and made of rubber material. Calca nea I epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf. Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.



**PH-8033**

### Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



## CORK INSOLE WITH ARCH SUPPORT 3/4 LENGTH

Longitudinal and metatarsal arch support manufactured from mixture of special cork and rubber material. Calcaneal epin and metatarsus head areas which take the toll of the body are softened. Lessens the pain at foot and calf, Arc support supplement for low base and related pain relief. Distributes the standing body weight evenly and reduces load and stress on the joints.



**PH-8034**

### Size:

S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



## STEEL INSOLE



Made of stout leather and sup-ported with stainless steel bridge has a longitudinal and metatarsal support. Distributes the standing body weight evenly and reduces load and stress on the joints. Rigid surface support does not allow the longitudinal arc to collapse.

**USAGE**  
Pes planus-planovalgus  
Pes cavus  
Traverse arc support

**PH-8035**

Size:	
S	35-37 cm
M	38-40 cm
L	41-43 cm
XL	43-45 cm



## SILICONE METATARSAL PAD WITH RING

It is made of fabric material that does not harm the skin. The position of the silicone ring is adjusted by inserting it on the 2nd finger. Metatarsals relieve pain by removing the pressure on their heads.

**USAGE**  
Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



**PH-8118**

**Size:**  
**Standard**

## 3/4 INSOLES (SILICONE)



Middle degree of hardness, is made from 100% silicon. 3/4 in size and offers comfortable use.

**USAGE**  
Provides anatomical support to the sole and enables even distribution of pressure on the sole. Has a relaxing and comforting effect.

**PH-8018**

Size:	
S	35-37 cm
M	38-40 cm
L	41-43 cm



## SILICON HEEL CUP

made of high quality of silicone material. Specially designed to fit figure of the foot. Softened calcaneal epin area lessens the load on the heels which causes pain. Lessens the pain at foot and calf.

**USAGE**  
Heel spur Plantar fasciitis  
Avascular necrosis of calcaneus (sever's disease)  
Heel pain



**PH-8019**

Size:	
S	34-37 cm
M	38-39 cm
L	40-41 cm



## SILICON TOE SEPERATOR



made of high quality of silicone material. Comfortable for daily use. Prevents deviation of the toe to the lateral in hallux valgus deformity.

**USAGE**  
Helps proper positioning of toes in the shoes and prevents overlapping.

**PH-8020**

**Size:**  
**Standard**

## SILICON WEDGE

It is made of silicone material and the upper surface is covered with fabric. In accordance with the problem, it is placed on the inner or outer side of the shoe. Due to its adhesive nature, it can be fixed within the shoe.

### USAGE

- 1-Calcaeus varus,
- 2-Calcaeus valgus,
- 3-0-Leg, X-Leg,
- 4-Meniscus degeneration mity.



**PH-8024**

### Size:

S	35-37 cm	
M	38-40 cm	
L	41-43 cm	
XL	43-45 cm	

## SILICONE HEEL CUP FABRIC COATED WITH EPINE



Manufactured from silicon material, covered with fabric. Softened calcaneal epin area prevents the load pressing onto a single point. Due to its adhesive nature, it can be fixed within the shoe.

### USAGE

- 1-Heel spur,
- 2-Plantar fasciitis,
- 3-Avascular necrosis of calcaneus (sever's disease),
- 4-Heel pain.

**PH-8025**

### Size:

S	35-37 cm	
M	38-40 cm	
L	41-43 cm	
XL	43-45 cm	

## SILICONE MATETARSAL PAD

It is made of silicone material and the upper surface is covered with fabric. It is placed within the proximal of the metatarsus heads. It relieves the pressure of the metatarsus heads and reduces pain. Due to its adhesive nature, it can be fixed within the shoe.

### USAGE

- 1-Over-sensitivity of metatarsus head,
- 2-Metatarsalgia,
- 3-Callus developing under the metatarsus,
- 4-Morton's foot-hallux rigi-dus-claw toe deformity,
- 5-Aceptycal necrosis of metatarsus head (Freiberg's disease).



**PH-8026**

### Size:

S	35-37 cm	
M	38-40 cm	
L	41-43 cm	
XL	43-45 cm	

## SILICONE ARC SUPPORT



It is made of silicone material and the upper surface is covered with fabric. It supports the longitudinal arch and reduces the increased tension in the sole. It starts from the calcaneus' inner tubercle and continues until the rear of the metatarsus heads. Due to its adhesive nature, it can be fixed within the shoe.

### USAGE

- 1-Pes planus,
- 2-Pes cavus,
- 3-When the longitudinal arch needs to be supported.

**PH-8027**

### Size:

S	34-37 cm	
M	38-39 cm	
L	40-41 cm	
XL	43-45 cm	

## SILICONE HAMMER TOE CUSHION (PAIR)

It is made of silicone material that does not harm the skin. It is used by passing on hammer fingers. It helps to eliminate deformation by applying support over the joint.

### USAGE

- Over-sensitivity of metatarsus head  
Metatarsalgia  
Callus developing under the metatarsus



**PH-8028**

**Standard Size:  
Left / Right**

## SILICONE METATARSAL PAD WITH RING

Manufactured of skin friendly, hypoallergenic silicon material, Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

**USAGE**  
Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



**PH-8016**

**Size:**  
**Standard**

## SILICONE METATARSAL PAD WITH TOESEPERATOR



Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

**USAGE**  
Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.

**PH-8029**

**Size:**  
**Standard**

## SILICON BUNION SHIELD

Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma, Used to reduce bunion sensitivity,

**USAGE**  
Manufactured of skin friendly, hypoallergenic silicon material. Prevents bunion friction and protects against trauma. Used to reduce bunion sensitivity.



**PH-8030**

**Size:**  
**Standard**

## SILICONE TOE CAP



Made of high quality of silicone material. Comfortable for daily use. Prevents deviation of the toe to the lateral in hallux valgus deformity.

**USAGE**  
Is made of silicone. Best to use in the daytime. Used to heal deviation in lateral toe Hallux Valgus.

**PH-8031**

**Size:**

S  
M





## VARICOSE SOCKS CALF (opened or closed toe)

### USAGE

Varicosis Stockings apply a decreasing compression towards the calf and upper leg whose external pressure is highest on the ankle.

Thus, they expedite the blood flow by providing a regular pressure distribution in legs and feet.

Patellar dislocations

Post-op recovery, protection and immobilization

### CLASS 1

Knitted elegantly and helps to relax  
Pressure offers of (14 mmHg)

**PH-8037**

### CLASS 2

Knitted elegantly and helps to relax  
Pressure offers of (20-30 mmHg)


**PH-8040**

### CLASS 3

Knitted elegantly and helps to relax  
Pressure offers of (30-40 mmHg)

**PH-8043**

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
cB	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
cB1	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
cC	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
cD	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
cE	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
cF	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
cG	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



Open Toe



## VARICOSE SOCKS THIGH (opened or closed toe)

### USAGE

Varicosis Stockings apply a decreasing compression towards the calf and upper leg whose external pressure is highest on the ankle.

Thus, they expedite the blood flow by providing a regular pressure distribution in legs and feet.

Patellar dislocations

Post-op recovery, protection and immobilization

### CLASS 1

Knitted elegantly and helps to relax  
Pressure offers of (14 mmHg)

**PH-8038**

### CLASS 2

Knitted elegantly and helps to relax  
Pressure offers of (20-30 mmHg)


**PH-8041**

### CLASS 3

Knitted elegantly and helps to relax  
Pressure offers of (30-40 mmHg)

**PH-8044**

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
<b>cB</b>	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
<b>cB1</b>	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
<b>cC</b>	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
<b>cD</b>	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
<b>cE</b>	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
<b>cF</b>	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
<b>cG</b>	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



Open Toe

## VARICOSE SOCKS (PANTY)

### USAGE

Varicosis Stockings apply a decreasing compression towards the calf and upper leg whose external pressure is highest on the ankle.

Thus, they expedite the blood flow by providing a regular pressure distribution in legs and feet.

Patellar dislocations

Post-op recovery, protection and immobilization

### CLASS 1

Knitted elegantly and helps to relax  
Pressure offers of (14 mmHg)

**PH-8039**

### CLASS 2

Knitted elegantly and helps to relax  
Pressure offers of (20-30 mmHg)


**PH-8042**

### CLASS 3

Knitted elegantly and helps to relax  
Pressure offers of (30-40 mmHg)

**PH-8045**

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
<b>cB</b>	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
<b>cB1</b>	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
<b>cC</b>	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
<b>cD</b>	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
<b>cE</b>	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
<b>cF</b>	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
<b>cG</b>	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



# RESİM YOK

## ANTIEMBOLISM SOCKS

### USAGE

prevent leg swelling and helps with leg discomfort, improves blood circulation, commonly known as blood clots, advised for post surgery or prolonged bed rest.

### ANTIEMBOLISM SOCKS CALF

Prevent varicose veins. It accelerates blood circulation. It increases the blood flow of the veins.



**PH-8046**

### ANTIEMBOLISM SOCKS THIGH

Prevent varicose veins. It accelerates blood circulation. It increases the blood flow of the veins.



**PH-8047**


### ANTIEMBOLISM SOCKS THIGH WITH BELT

Muscles pain. Neck arthritis. Sliding discs. After the traumas abd bruises.



**PH-8048**

Leg Circumference (cm)	0	1	2	3	4	5	6	7	8
<b>cB</b>	16-18	18-20	20-22	22-24	24-26	26-28	28-30	30-32	32-35
<b>cB1</b>	21-25	23-27	24-29	26-32	29-35	31-37	33-39	35-41	37-44
<b>cC</b>	26-31	28-34	30-37	33-40	35-43	37-45	39-47	41-49	43-52
<b>cD</b>	25-30	27-33	29-36	32-39	34-42	36-44	38-46	40-48	42-51
<b>cE</b>	29-35	32-38	34-41	37-44	39-47	41-49	43-51	45-53	47-55
<b>cF</b>	34-43	37-47	40-52	43-56	46-60	50-65	53-68	56-71	60-78
<b>cG</b>	40-53	43-57	46-62	49-67	52-72	56-76	60-79	64-82	69-92



# COMPRESSION STOCKINGS



# PEDIATRIC PRODUCTS



## SOFT CERVICAL COLLAR FOR KIDS

It provides superior comfort and support to the neck through anatomical designed. It refrains from excessive cervical activity and provides extra warmth around the neck. Compatibility MR.

### USAGE

Muscles pain. Neck arthritis. Sliding discs. After the traumas and bruises



**PH-9001**

**Size:  
Standard**

## CERVICAL COLLAR (NELSON) FOR KIDS

RESİM YOK

It provides superior comfort and support to the neck through anatomical designed. It helps the patient during the treatment by restricting the movements in the neck area. Compatibility MR.

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
After the traumas and bruises

**PH-9002**

**Size:  
Standard**

## CERVICAL COLLAR (PHILEDELPHIA) FOR KIDS

It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It has a chin support. Compatibility MR.

### USAGE

Soft tissue trauma.  
Ligaments and bones injuries.  
gliding discs.  
Stabilizations in the cases of the unconsciousness.

RESİM YOK

**PH-9003**

**Size:  
Standard**

## CERVICAL COLLAR (PHILEDELPHIA-TRACHEOSTOMY) FOR KIDS

RESİM YOK

It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It is a collar with trachea open. Compatibility MR.

### USAGE

Soft tissue trauma.  
Ligaments and bones injuries.  
Sliding discs.  
Stabilizations in the cases of the unconsciousness.

**PH-9004**

**Size:  
Standard**

## CERVICAL COLLAR WITH CHIN SUPPORT (PLASTIC) FOR KIDS

Made of PVC with chin support in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine.

### USAGE

Muscles pain.  
Neck arthritis.  
Sliding discs.  
Cases of radiculopathy.  
After the traumas and bruises.



**PH-9005**

**Size:**  
**Standard**

## CERVICAL COLLAR (PLASTIC) FOR KIDS



Made of PVC in a design that can be modified as wanted to stabilize the neck and reduce the pressure on the spine

### USAGE

Sliding discs  
Neck arthritis  
After the traumas and bruises

**PH-9006**

**Size:**  
**Standard**

## ARM SLING MESH FABRIC FOR KIDS

manufactured from perforated fabric. Products design provides patient's arm to be fixed parallel to the near trunk. It can be used for both sides.

### USAGE


Fractures and bruises of the arm.  
Bruised shoulder.



**PH-9007**

### Size:

S	16-22 cm
M	22-28 cm
L	28-34 cm



## CLAVICULA BANDAGE FOR KIDS



Made of comfortable cloth covered from the inside with a layer of soft cloth which does not cause allergies, helps to stabilize the wishbone thanks to belts that pull back the shoulders to the extent required to give the appropriate positioning.

### USAGE

Clavicle dislocations.  
Clavicle fractures.

**PH-9008**

### Size:

S	45-50 cm
M	50-60 cm
L	60-70 cm





## ELBOW SUPPORT FOR KIDS

Made of neoprene fabric that provides warmth and protection. Thanks to the three-headed silicone pad inside, it applies constant pressure and compression to the epycondil heads.

**USAGE**  
Lateral epicondylitis.  
Internal articular fluid.



**PH-9009**

Size:

S	16-22 cm
M	22-28 cm
L	28-34 cm



## WRIST BANDAGE FOR KIDS



Provides protection for the wrist, sprains and bruises wrist  
Retains body heat to improve bloody circulation and accelerate wound healing.

**USAGE**  
Twists, Trumas, bruises.  
Protection while doing sports,  
Protection at work when the jobs involves using the wrist a lot.

**PH-9010**

Size:  
**Standard**

## WRIST SPLINT FOR KIDS

Made of neoprene fabric that provides warmth and protection. Constrains movements of the wrist thanks to the aluminium support inside.

**USAGE**  
Wrist tunnel cases.  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.



**PH-9011**

Size:

S	9-11 cm
M	11-13 cm
L	13-15 cm



## WRIST SPLINT WITH THUMB FOR KIDS



There is an aluminium support under the wrist and near the thumb. It can use to both hands by changing aluminium support under the wrist. There is a velcro that increases the stabilization around the wrist.

**USAGE**  
Wrist tunnel cases,  
Tendonitis, tenosynovitis.  
After taking off gypsum.  
After Surgeries.

**PH-9012**

Size:

S	9-11 cm
M	11-13 cm
L	13-15 cm





## NEOPREN KNEE SUPPORT WITH PAD FOR KIDS

Made of neoprene fabric that provides warmth and protection. Used during doing sports.

### USAGE

Supporting the Knee and protecting it from injuries.



**PH-9013**

### Size:

S	25-28 cm
M	28-31 cm
L	31-34 cm



## NEOPREN KNEE SUPPORT WITH LIGAMENTS FOR KIDS



Made of neoprene fabric that provides warmth and protection. The flexible underNire on the sides increase stabilization and thus decrease the load on the knee. The soft neoprene blend retains body heat to improve circulation and promote healing.

### USAGE

Arthritis, bursitis. Increase of intra articular liquid. Pain induced by knee cap  
Meniscus problems. 1ST degree twist and ligament instability, Post-op or trauma for

**PH-9014**

### Size:

S	25-28 cm
M	28-31 cm
L	31-34 cm



## PEDIATRIC KNEE IMMOBILIZER FOR KIDS

It provides immobilization in extension. It is designed anatomically and includes aluminium supports to provides stabilization in the posterior, medial and lateral. Thanks to the three-piece pad design, it can be easily adjusted for leg circumference.

### USAGE

Recovery after ligament injury. Meniscus Post-op recovery. Knee joint in trauma cases Ostrochondra I fractures. Acute prainsstrains of the knee. Patellar dislocations. Post-op recovery, protection and immobilization



**PH-9015**

### Size:

S	25-28 cm
M	28-31 cm
L	31-34 cm



## ANKLE SUPPORT FOR KIDS



Made of neoprene fabric that provides warmth and protection.. Velcro fasteners around the ankle provide ease at wear, Supports the soft neoprene blend retains body heat to improve circulation and promote healing process.

### USAGE

1st degree twist and ton us  
Minor fractures  
Compression on edema  
Protection during and after the sportive activities

**PH-9016**

### Size:

S	18-20 cm
M	20-22 cm
L	22-24 cm



## FIGURE OF EIGHT BANDAGE FOR KIDS

Made of neoprene fabric  
Contains pores that help for ventilation. The Crossed lacing provides extra pressure and stabilization thanks to the rubber pads on the sides.

### USAGE

Ligament instability and rupture  
Edema prevention  
Post-op or trauma for protection and rehabilitation purposes



**PH-9017**

**Size:**  
**Standard**

## PAVLIK SUPPORT FOR KIDS



Made of 3 cm cotton belts, non-elastic, stabilize the hip and the feet in a perfect position.

### USAGE

Congenital Dislocation

**PH-9018**

### Size:

S	0-3 Month
M	3-6 Month
L	6-9 Month

## FREJKA ORTHOSIS FOR KIDS

Manufactured from 12 mm top quality molded plastazote foam.  
Holes allow air to circulate, soothing support and relief.

Latex-free, non-toxic and hypoallergenic plastizode foam material will not cause skin irritation. Velcro straps are used for easy donning and doffing. Plastizode foam material also is radiolucent which allows X-ray, it possible to use CT, MRI. Positions at the hip joint in flexion and abduction.

### USAGE

Hip subluxations



**PH-9019**

### Size:

S	0-3 Month
M	3-6 Month
L	6-9 Month

## IMBILICAL HERNIA TRUSS FOR KIDS



Made of flexible cotton fabric. The pad applies pressure to the inguinal hernia. Perineal strap prevents sliding. It comes in standard size but adjustable.

### USAGE

Umbilical hernia

**PH-9020**

**Size:**  
**Standard**

### INGUINAL RAPTURE BELT FOR KIDS

Made of a special cotton fabric. Does not effect the daily life. Provides support after surgeries. Fits perfectly under clothes.

**USAGE**  
Single sided inguinal hernia



RESİM YOK

**PH-9021**

**Size:**  
Standard

### INGUINAL RAPTURE PEDIATRIC BELT (DOUBLE) FOR KIDS



Made of a special cotton fabric. Does not effect the daily life. Provides support after surgeries. Fits perfectly under clothes.

**USAGE**  
Double sided inguinal hernia

**PH-9022**

**Size:**  
Standard

### DENNIS BROWN (MOVABLE/STABLE) FOR KIDS

The child can do independent movement of the legs such as kicking or crawling. decrease the risk of the boot sliding off of the foot and regional rash. The child will be active and sleeps more comfortably.


**USAGE**  
Used for hip dislocation and pes equinovarus deformity.



**PH-9023**

**Size:**  
Standard

### DENNIS BROWN SHOES FOR KIDS



RESİM YOK

**PH-9024**

**Size:**  
From 17  
to 27



# OTHER PRODUCTS





## SILICON BREAST

It is used to provide a natural appearance when chest tissue is removed after mastectomy, to reduce pain caused by muscle imbalances, to protect the surgical wound from external bumps and crushes. The laundry is made of cotton and non-allergenic special fabric and silicone material. Shoulder straps are softened for comfortable use. The pockets on both sides allow the prosthesis to be used as right and left.



**RESİM YOK**

**PH-100**

**Size:  
Standard**

## SILICON BREAST'A BRA



**RESİM YOK**

The prosthetic bra, which is very comfortable and stylish to wear, is specially designed for breast prosthesis. It supports the prosthesis and weight and provides a sense of confidence and comfortable use. It is made of cotton fabric, the inner surface of which does not irritate the skin. Shoulder straps are softened for comfortable use. The pockets on both sides allow the prosthesis to be used both left and right.

**PH-101**

**Size:  
Standard**

## PORTABLE CLOSED SEAT LIFTER (WITHOUT COVER)

Made of 3 cm cotton belts, non-elastic, stabilize the hip and the feet in a perfect position.

**USAGE**  
Congenital Dislocation



**RESİM YOK**

**PH-102**

**Size:  
Standard**

## PORTABLE CLOSED SEAT LIFTER (WITH COVER)



**RESİM YOK**

Latex-free, non-toxic and hypoallergenic plastizode foam material will not cause skin irritation. Velcro straps are used for easy donning and doffing. Plastizode foam material also is radiolucent which allows X-ray, it possible to use CT, MRI. Positions at the hip joint in flexion and abduction.

**USAGE**  
Hip subluxations


**PH-103**

**Size:  
Standard**

## SILICON EXERCISE BALLS

Made of high quality silicone material which does not cause allergies and has no smell. Designed to withstand the pressure applied on, without losing their shape. Available in three different colors.

**USAGE**  
Used so as to rehabilitation. Three different hardnesses existing. Red - Hard, Blue - Medium, Green - Soft



**RESİM YOK**

**PH-104**

**Size:  
Standard**

VISCO NECK PILLOW

RESİM YOK

Size:  
60\*40\*11/8 cm

PH-200

VISCO NECK PILLOW

RESİM YOK

Size:  
60\*37\*11/9 cm

PH-201

VISCO CLASSIC  
NECK PILLOW

RESİM YOK

Size:  
60\*40\*12 cm

PH-202

VISCO CLASSIC  
NECK PILLOW

RESİM YOK

Size:  
60\*40\*15 cm

PH-203

SNORING REDUCING  
VISCO PILLOW

RESİM YOK

Size:  
55\*40\*12/10 cm

PH-204

APNEA PILLOW

RESİM YOK

Size:  
55\*40\*12/10 cm

PH-205

VISCO NECK  
PILLOW WITH GEL

RESİM YOK

Size:  
57\*37\*11/7 cm

PH-206

VISCO CLASSIC NECK  
PILLOW WITH GEL

RESİM YOK

Size:  
60\*40\*15 cm

PH-207

# PILLOWS

VISCO NECK  
PILLOW WITH GEL

RESİM YOK

Size:  
60\*37\*11/9 cm

PH-208

FLAT CLASSIC  
PILLOW WITH GEL

RESİM YOK

Size:  
55\*40\*11 cm

PH-209

VISCO TRAVEL  
PILLOW

RESİM YOK

Size:  
31\*27\*7 cm

PH-210

VISCO TRAVEL  
PILLOW

RESİM YOK

Size:  
31\*27\*10 cm

PH-211

VISCO SEAT  
CUSHION

RESİM YOK

Size:  
43\*9 cm

PH-212

VISCO SEAT  
CUSHION (OPEN)

RESİM YOK

Size:  
43\*43\*9,5 cm

PH-213

VISCO SQUARE  
CUSHION

RESİM YOK

Size:  
45\*43\*7 cm

PH-214

VISCO U SQUARE  
CUSHION

RESİM YOK

Size:  
41\*39\*6 cm

PH-215

VISCO ANATOMIC SEAT

RESİM YOK

Size:  
42\*40\*9/2 cm

PH-216

VISCO COCCYX CUSHION

RESİM YOK

Size:  
44\*36\*9,5 cm

PH-217

VISCO POSTURE SEATING CUSHION

RESİM YOK

Size:  
44\*38\*13/7 cm

PH-218

VISCO SQUARE CUSHION WITH GEL

RESİM YOK

Size:  
45\*43\*7 cm

PH-219

VISCO ANATOMIC SEAT WITH GEL

RESİM YOK

Size:  
42\*40\*9/2 cm

PH-220

BETWEEN LEG VISCO PILLOW

RESİM YOK

Size:  
25\*19\*14 cm

PH-221

BACK AND WAIST SUPPORT PILLOW

RESİM YOK

Size:  
54\*42\*9 cm

PH-222

VISCO LUMBAR CUSHION

RESİM YOK

Size:  
38\*32\*12 cm

PH-223



# PILLOWS

VISCO LUMBAR  
CUSHION

RESİM YOK

Size:  
35\*34\*12 cm

PH-224

VISCO LUMBAR  
CUSHION

RESİM YOK

Size:  
34\*32\*8 cm

PH-225

VISCO LUMBAR  
CUSHION

RESİM YOK

Size:  
39\*34\*9 cm

PH-226

VISCO REFLUX  
PILLOW

RESİM YOK

Size:  
65\*60\*24 cm

PH-227



**PhysiPlus**<sup>®</sup>  
REHAB & SPORT

**İSTANBUL ŞEHİR MEDİKAL DEPOSU İMALAT VE DIŞ TİC. LTD ŞTİ.**

Kayabaşı Mah. Ulubatlı Hasan Cad.

Kayaşehir C10 Blok No:2L/2 Başakşehir / İstanbul

T: +90 212 809 80 80 - +90 541 534 89 34

[www.istanbulsehirmedikal.com](http://www.istanbulsehirmedikal.com)