







NINGBO TECHART MEDICAL EQUIPMENT CO.,LTD

Add: No.99 South Hongtang Road Ningbo China 315033

Tel: +86-574-83076915 Fax: +86-574-87137656

E-mail:sales01@techartmed.com

sales02@techartmed.com sales03@techartmed.com

sales03@techartmed. Web: www.techartmed.com Your Ideal Operating Table Supplier

> Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery






Modular Design

The modular structure of the table top is tailored for various surgical discipline needs ,which also can low down the cost of the hospital.

Easy Click Connection Design

The easy click connection design offers a convenient, safe, intuitive modular mounting solution. Different surgery have different requirement of table. The easy click connection make the table easy to adapt to the surgeries.

Two Joint Head Plate



is controlled by manual handle.

TECHARTMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery





TECHARTMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery

06

Multi-version tabletop -Flexible configuration

In addition to the standard configuration of the table, ET800 also offers a variety of table combinations for customers to choose from flexibly, to meet the needs of different procedures





ET 800 be extended and configured with different surgical accessories

ET800 with standard tabletop



Back plate in 2 sectios



Double control leg plate Segmental adjustment of thigh and lower leg



Tabletop Version 1





Motorised and separated leg plate

Back plate in two sections



Tabletop Version 2






Separated manual E control leg plate

Back plate in two sections



Tabletop Version 3








Back plate in one section



Tabletop Version 4







Back plate in one section



TECHARTMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery





Extreme High Extreme Low Position

Maximum and minimum heights of the ET800, not only suitable for kinds of an surgery also for doctors' comfortable working, both seating and standing.



TECHARTMED®

Flexible Patient Position

ET800 is a multi-functional operating table that meets various requirements, such as chest, kidney, cholecyst surgery, Endoscopic surgery, gynecologic surgery, ENT SURGERY, etc.

08



Optional Position Mattress, Leg Support, Waste Basin



General Surgery



Lateral Lying Position For Chest Kidney, cholecyst Operation



E.N.T. Operation Position



Chest And Abdomen Operation Position



Optional Lateral Arm Support



Anorectum Operation Position

Model Site	length	Width	Tabletop lowest height	Elevation range
ET800	≥2100mm	≥550mm	550mm	≥500mm

Ps1: without side rail.

Ps2: lowest height without mattress and x-ray board.

PARAMETER ET800 Electro-Hydraulic Operating Table

Movements		Model:ET800	
Tobleton transcionh::	Reverse Trendelenburg	≥35°	
Tabletop trendelenburg	Trendelenburg	≥35°	
Lateral Tilt (Left And Right)	≥25°	
	Up	≥50°	
Head plate	Down	≥60°	
	Flexible	≥50°	
Back plate	Up	≥80°	
Back plate	Down	≥40°	
	Electric Up	≥80°	
Leg plate	Electric Down	≥90°	
Leg plate	Manual Down	≥90°	
	Manual Outward	≥90°	
Kidney Brid	lge Elevation	120 mm	
Horizon	350 mm		
FI	≥220°		
Re	≤100°		


OPTIONAL UNIT

Hand Remote	Standard	User-defined memory location function	Standard	Extension modules	Standard (1 piece) Optional more
Control Panel	Standard	Extreme high extreme low position	Standard	Electric leg plate (one piece)	Optional
Wireless LCD hand remote	Optional	Kidney bridge	Standard	Electric split leg (independent Controlled)	Standard


TECHARTMED

Your Ideal Operating Table Supplier


Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery



Horizontal Sliding:350mm



Kidney Bridge Elevation:120mm








Table Highest Height: 1050mm













Head plate flex

Head plate up

Head plate down














Lateral tilt


Trendelenburg position














Leg plate horizontal position