



MY·FILLER
HA DERMAL FILLERS





New Life to the Skin

TABLE OF CONTENTS:

4	Skinboosters
6	HA Refill
8	HA Retouch
10	Dual-Link Technology
12	Glips
14	Max
16	Trio Technology
18	Soft
20	Medium
22	Strong
24	Range Comparison
26	My Med Aesthetics
28	Beauty Journey

New Life to the Skin

WITH OUR HA DERMAL FILLERS

My Filler features:

- Two different technologies: Dual-Link and TRIO.
- Unique and innovative solutions with 7 different formulations.
- EU certification by an internationally recognized Italian Notified Body.
- Several MOH International registrations.
- Ergonomically designed finger grip and comfortable plunger.
- Beauty Pass inside for your patients' convenience.
- Maximum quality, safety and efficacy in the clinical practice.
- Designed to enhance natural beauty.
- Dynamic formulation with perfect balance between elasticity and viscosity.







THE NEXT GENERATION **Skinboosters**

Standard skinboosters are simply non-cross linked HA, which hydrates the skin with short action time.

Our dynamic 2nd generation skinbooster is a real pioneer in bringing the concept of skinboosters to the next level.

The secret is a big touch of cross-linked hyaluronic acid for the long-lasting result combined with linear hyaluronic acid for the immediate effect to increase the moisturizing factor even further.







NIX FILLER
100 ml



HA REFILL

Enhanced to a New Level

Our pioneer 2nd generation skinbooster combines cross-linked HA and linear HA for super moisturized skin in the long run.

- **Recommended areas:** the entire face.
- **The result:** a hydrated, radiant-looking and even-textured complexion, that lasts longer with increased effectiveness.
- **Other characteristics:** low viscosity with the dynamic material.
- **Extra characteristic:** excellent for décolleté and neck rejuvenation as well.

HA Concentration	Cross-linked rate	Tan δ
15+5 mg/ml	1,5	0,45



MY-FILLER
HA RETOUCH

G

HA RETOUCH Boost Even Further

Double dose of our unique 2nd generation skinbooster for an immediate and long-lasting result.

- **Recommended areas:** the entire face.
- **The result:** a hydrated, radiant-looking and even-textured complexion, that lasts longer with increased effectiveness and a double dose.
- **Other characteristics:** low viscosity in a syringe of 2 ml.
- **Extra characteristic:** excellent for décolleté and neck rejuvenation as well.

HA Concentration	Cross-linked rate	Tan δ
15+5 mg/ml	1,5	0,45



IMMEDIATE AND LONG-LASTING **Dual-Link Technology**

The Dual-Link Technology was created as a result of years of development and consultations with specialists and patients, which provides a solution to their wishes: the achieved result is not only long-lasting but can be seen and felt immediately!

The Dual-Link Technology combines the advantages of both linear HA and cross-linked HA in one unique line of products. The special formulation provides an immediate hydration effect, thanks to the linear HA, as well as a long-lasting result due to the cross-linking HA, ensuring patient satisfaction even further.

The crossed-linked and linear HA combination makes the formula more dynamic because it gives relatively high elasticity even with a lower viscosity. Therefore it can be easily shaped for superficial treatments (My Filler Glips) and easier to inject for deep treatments (My Filler Max) as well.







MY-FILLER
GLIPS



GLIPS Long-Lasting Glowing Lips

My Filler Glips is dedicated and exclusively created for lip treatments. Enhance lips with hydration, added volume and re-contouring feature.

- **Recommended injection area:** lips.
- **The result:** immediate, long-lasting, shining lips with natural-looking .
- **Other characteristics:** medium viscosity and dynamic formulation.

HA Concentration	Cross-linked rate	Tan δ
20+5 mg/ml	2	0,27



MY-FILLER
MAX




MAX Unique Solution Maximized Volume

My Filler Max gives maximum volume with the immediate shining of the skin. With its high HA concentration can be an excellent choice for deep-volume correction.

- **Recommended injection area:** cheeks, chin, jaw, nose and midface area.
- **For men:** a more masculine appearance can be achieved with a strong and well-defined jawline, shape and volume.
- **For women:** can be good for a strong volume correction with immediate shining effect.
- **Other characteristics:** high viscosity with excellent formability.

HA Concentration	Cross-linked rate	Tan δ
22+8 mg/ml	3	0,18



SCIENTIFIC INNOVATION **Trio Technology**

Three dermal fillers are based on TRIO Technology. Soft, Medium and Strong formulas contain specific concentrations of hyaluronic acid combined with different cross-linking rates in a special manufacturing environment. As a result, three fillers with largely different properties are obtained allowing a wide range of usability.

WHY DO PROFESSIONALS AND PATIENTS CHOOSE OUR 'BIG THREE'?

- High biocompatibility.
- Personalised solutions tailored to your patient's needs.
- High hydrophilic capacity.
- Easy to shape.
- Long-lasting effect.
- Stability after injection.
- Innovative solutions for various treatment.
- Safe to inject, with no skin reaction post-treatment.







SOFT Fine Lines General Hydration

My Filler Soft is recommended for situations where the area to be treated requires the elimination of fine wrinkles with a slight increase of volume and super hydration.

Excellent choice of prevention and for delaying the formation of wrinkles, and for a hydrated healthier appearance.

- **Recommended areas:** Periorbital (tear troughs), perioral and lips.
- **Other characteristics:** low viscosity and superficial injection possibilities enable to enhance the results of various treatments with its long-lasting moisturizing effect.

HA Concentration	Cross-linked rate	Tan δ
20 mg/ml	1,5	0,46





MY FILTER
MEDIUM

ARABIAN

MEDIUM

Moderate Volume

Energetic Hydration

Thanks to its universal application properties, My Filler Medium has become one of the most popular products from the My Filler product portfolio. It can be used perfectly if our goal is to achieve medium volume and energetic hydration in the treated area.

- **Recommended areas:** nasolabial wrinkles, mentolabial (marionette wrinkles), glabella and lip volume increase.
- **Other characteristics:** the viscosity of the material is suitable for volume correction and for re-shaping as well.

HA Concentration	Cross-linked rate	Tan δ
25 mg/ml	2,5	0,30





MY FILTER
STRONG

STRONG Innovative High-Volume

My Filler Strong's high HA concentration makes it outstanding in the world of hyaluronic acid fillers. Recommended especially for contouring and high-volume corrections.

- **Recommended areas:** Cheek, chin, jawline, nose and midface.
- **Other characteristics:** high viscosity as well as excellent formability properties. After injection the material remains stable in the area avoiding any sort of unexpected migration.
- **Extra characteristic:** can be well combined with My Filler Medium or even My Filler Max with the well-known 'sandwich technique' for a multilayered approach that comprehensively treats the volume loss.

HA Concentration	Cross-linked rate	Tan δ
30 mg/ml	3,5	0,16



Range comparison

	SKINBOOSTERS		
	HA REFILL	HA RETOUCH	GLIPS
HA TYPE	Monophasic	Monophasic	Monophasic
CROSS-LINK AGENT	BDDE	BDDE	BDDE
CROSS-LINK RATE	low	low	moderate
HA CONCENTRATION	15+5 mg/ml	15+5 mg/ml	20+5 mg/ml
SYRINGE VOLUME	1 ml	2 ml	1 ml
SUGGESTED NEEDLE	27G, 30G, 32G	27G, 30G, 32G	27G, 30G, Cannula
PRODUCT CODE	MYF107	MYF108	MYF105
INJECTION DEPTH	Epidermis	Epidermis	Derma



DERMAL FILLERS			
MAX	SOFT	MEDIUM	STRONG
Monophasic	Monophasic	Monophasic	Monophasic
BDDE	BDDE	BDDE	BDDE
high	low	moderate	high
22+8 mg/ml	20 mg/ml	25 mg/ml	30 mg/ml
1 ml	1 ml	1 ml	1 ml
25G - 27G, Cannula	27G, 30G, 32G	27G, 30G, Cannula	25G - 27G, Cannula
MYF104	MYF001	MYF002	MYF003
Subcutaneous	Epidermis	Derma	Subcutaneous



MY·FILLER
HA DERMAL FILLERS



A PIONEER WITH GLOBAL MISSION.
Unique solutions.

My Med Aesthetics is the manufacturer of three premium brands: **My Filler** - professional aesthetic dermal fillers, **Revitalize** - skincare professional serums and **Oh!** - skincare facial creams. Our mission is to spread the message of conscious slow-aging in harmony and synergy between body and soul.

Founded in 2016 and nowadays represented in more than 50 countries worldwide, **My Med Aesthetics** is empowering medical aesthetics practitioners and their patients with clinically proven, high-end innovative solutions for skin rejuvenation.

Certified company by an international reliable Notify Body, **My Med Aesthetics** obtained the UNI EN ISO 9001:2015 for Quality Management System and the UNI EN ISO 13485:2016 for Manufacturing System. Moreover, our brand received the Certificates of Registration from the EUIPO (European Union Intellectual Property Office), protecting our partners and guaranteeing their work safety. Moreover, the product certification is done by one of the most important Notify Body from Italy, Istituto Superiore di Sanita, constantly auditing our company supervising the respect of all rules and high standards of quality.



DEDICATED TO ELEVATE YOU AND YOUR PRACTICE!

Rethinking medical aesthetics as a holistic approach for inner and outer beauty we are devoted to constantly grow together with our aesthetic partners to enhance natural beauty and find harmony in each individuals' **Beauty Journey**.


Our global presence



Evoke your Beauty Journey!

Have you realized that you have your own Beauty Journey? Which has its own ups and downs but leads you to the direction of acceptance and perfect balance?

My Med Aesthetics supports you in this eternal cycle in order to restore, reboot and reshape the beauty of your skin and reconnect to your true self even more.



WWW.BEAUTYJOURNEY.EU

Oh!

Oh! - Make magic happen. A professional skincare range with simplicity and state-of-the-art technology at its core. The range contains key essentials medically proven to support all skin types at any age, boosted by My Med Aesthetics' unique and innovative **Epidermal Growth Factor**. With regular use, your clients can nourish their own skin at home through an easy-to-follow daily skincare routine using effective ingredients. Improve hydration within the layers of the skin, delay the signs of aging and help them optimize the results they have had through professional treatments.

RE
VITA
LIZE

Revitalize - Realize your inner beauty. The **Revitalize** 100% organic range gives the practitioner the ability to address multiple issues within the skin or the face, body, and scalp. Improve the texture and tone of the skin offering deep hydration, nourishment and a radiant complexion through the use of proven vitamins, amino acids and minerals. Reduce the appearance of fine lines and wrinkles while stimulating the production of collagen. Make skin irregularities disappear by promoting the growth of new skin cells with our extraordinary **Epidermal Growth Factor** as well as treating acne and scarring. Promote hair growth by increasing blood flow to the scalp and stimulating hair follicles the list goes on. Create cocktails with active ingredients to get the most out of the benefits that the **Revitalize** range can offer to your clients. By improving overall skin health and appearance, your clients will experience a boost to their confidence and self-esteem.

MY·FILLER

My Filler - New life to the skin. The cherry on the top when you need to add volume, smooth out wrinkles, contour or enhance facial features on a professional level. **My Med Aesthetics'** award-winning dermal filler line with its unique **Dual-link technology** contains monophasic cross-linked and non-cross-linked hyaluronic acid with second-generation skin boosters for a natural and long-lasting outcome.

Reconnect to natural beauty! Healing ourselves on a deeper level elevates the outer results. Essential to always remember: as in, so outside! The Beauty Journey encourages clients to maximize their results by combining the full range of My Med Aesthetics products for optimum outcome and customer satisfaction.

Version 2.0/2024



www.myfiller.eu



[mymedaesthetic](https://www.facebook.com/mymedaesthetic)



[@mymedaesthetics](https://www.instagram.com/mymedaesthetics)



[MyMedVideos](https://www.youtube.com/MyMedVideos)



Téglagyári út 5-11. 2097 Pilisborosjenő, Hungary
info@my-med.eu - www.my-med.eu

