



TRD REVISION KNEE SYSTEM



www.tipsan.com.tr

*the key in
orthopaedic
sciences*

TRD REVISION FEMORAL COMPONENT POSTERIOR STABILIZED (PS)

- CoCrMo Alloyed
- 6 Size, Left and Right
- Cemented



TRD REVISION EXTENSION STEM

- Titanium Alloyed
- To be Used for Both Femur and Tibia
- Diameters :
11 mm - 13 mm - 15 mm
17 mm - 19 mm - 21 mm
- Length :
65 mm
85 mm
105 mm
125 mm
150 mm



TRD REVISION FIXED TIBIAL COMPONENT

- CoCrMo Alloyed
- 6 Size Options
- Cemented



TRD REVISION OFFSET STEM ADAPTER


- Titanium Alloyed
- Offset : 0 - 3 - 5 mm



Fixation Screw to Offset Stem Adapter

TRD REVISION FIXED INSERT POSTERIOR STABILIZED (PS)

- UHMWPE High Crosslink
- 6 Size Options
- 9 Different Thickness Size are Available
8 mm - 10 mm - 12 mm
14 mm - 17 mm - 20 mm
23 mm - 26 mm - 29 mm



Fixation Screw to Tibia Component



TRD REVISION PATELLA COMPONENT

- UHMWPE
- 6 Size
- Diameters :
5 mm - 10 mm



TRD REVISION FEMORAL POSTERIOR BLOCK AUGMENT

- Titanium Alloyed
- 6 Size
- Thickness :
5 - 10 mm



Fixation Screw to Femoral Component



TRD REVISION FEMORAL DISTAL BLOCK AUGMENT

- Titanium Alloyed
- 6 Size
- Thickness :
4 - 8 - 12 mm



Fixation Screw to Femoral Component



TRD REVISION TIBIAL BLOCK AUGMENT

- Titanium Alloyed
- Thickness :
5 - 10 mm
- 6 Size Options



Fixation Screw to Tibial Metal Back

