

DERMAL FILLERS

PRODUCT RANGE OVERVIEW

Hafiller Model	Fine Lines	Derm	Derm Deep	Derm Plus	Sub Skin
Indication	Fine Lines around eye and mouth. Such as crow's feet and perioral lines.	Moderate wrinkles and skin folds. Such as frown lines. Lip enhancement.	Deep Wrinkles and folds. Such as nasolabial fold. Lip enhancement and face filling.	Breast and buttock augmentation.	Correction of facial contours include cheek, chin and nose.
Injection site					
Particle size	0.10-0.15mm	0.15-0.28mm	0.28-0.50mm	0.50-1.25mm	1.25-2.0mm
Concentration	20mg/ml	20mg/ml	20mg/ml	20mg/ml	20mg/ml
Needle size	30G	30G	27G or 26G	23G or 25G	23G or 25G
Duration	6-9 months	6-12 months	6-12 months	9-18 months	9-18 months