
# INTRODUCER, BOUGIE, GUIDE, EXCHANGE 插管导入器




#### **FEATURES AND BENEFIT:**

- The Introducer / Bougies used to facilitate endotracheal intubation.
- Made of medical grade PE, LATEX Free, Flexible.
- Curved tip can be passed blindly into the trachea when visualization of the glottis is inadequate to guide placement of endotracheal tube.
- Distal rounded smooth tip allows the bougies to be inserted without damaging the soft tissue.
- Memory retention allows shaping to individual preferences and needs.
- Clearly graduated marking provides depth indication to ensure correct placement.
- Have different sizes, suitable for pediatric, child and adult.





• How to know if intubation is expected to be difficult or not ?



Bougies device used in case of difficult intubation (Cormack superior to 2):

- Standard procedure recommendation for difficult intubation is to use as first intent long and flexible guides.
- This technique is very effective. According to several studies, the bougies treat more than 96% of difficult intubations.
- The bougie is the first intent device in the case of difficult intubation.

We have many ranges of Endotracheal tube Bougies.

## 1. Tracheal Tube Introducer (Bougies), Curved tip

Curved tip

		<b>1</b> 10	•25	■30	■35	■ 40
REF number	Size	OD	Length	Recommended ETT size	_	
PA210105	5Fr	1.7mm	53.5cm	ID 2.0mm - 3.0mm		
PA210106	6Fr	2.0mm	53.5cm	ID 2.5mm - 3.5mm		
PA210110	10Fr	3.3mm	70cm	ID 4.0mm - 5.5mm		
PA210114	14Fr	4.7mm	70cm	ID 6.0mm - 10.0mm	_	
PA210115	15Fr	5.0mm	70cm	ID 6.0mm - 10.0mm		
					_	

<sup>\* 10</sup>pcs/box, 200pcs/ctn

# 2. Tracheal Tube Introducer (Bougies), Straight tip

Straight tip

		■ 10	●25	■ 30	■ 35	■ 40
REF number	Size	OD	Length	Recommended ETT size	_	
PA210205	5Fr	1.7mm	53.5cm	ID 2.0mm - 3.0mm		
PA210206	6Fr	2.0mm	53.5cm	ID 2.5mm - 3.5mm	_	
PA210210	10Fr	3.3mm	70cm	ID 4.0mm - 5.5mm		
PA210214	14Fr	4.7mm	70cm	ID 6.0mm - 10.0mm	_	
PA210215	15Fr	5.0mm	70cm	ID 6.0mm - 10.0mm		

<sup>\* 10</sup>pcs/box, 200pcs/ctn

# 3. Tracheal Tube Introducer (Bougies), Easy-Carry Pack

Easy-Carry

FEATURES: Innovative compact packaging designed specifically for use by EMS professionals.




Tip	REF number	Size	OD	Length	Recommended ETT size
Curved tip	PA210510	10Fr	3.3mm	70cm	ID 4.0mm - 5.5mm
Curved tip	PA210514	14Fr	4.7mm	70cm	ID 6.0mm - 10.0mm
Curved tip	PA210515	15Fr	5.0mm	70cm	ID 6.0mm - 10.0mm
Straight tip	PA210510S	10Fr	3.3mm	70cm	ID 4.0mm - 5.5mm
Straight tip	PA210514S	14Fr	4.7mm	70cm	ID 6.0mm - 10.0mm
Straight tip	PA210515S	15Fr	5.0mm	70cm	ID 6.0mm - 10.0mm

<sup>\* 10</sup>pcs/box, 200pcs/ctn

## 4. Tracheal Tube Introducer (Bougies), Vented tip, without connector


Vented tip



<sup>\* 10</sup>pcs/box, 200pcs/ctn

## 5. Tracheal Tube Introducer (Bougies), Vented tip, with TPE connector


Vented tip



<sup>\* 5</sup>pcs/box, 100pcs/ctn

# 6. Tracheal Tube Introducer (Bougies), Vented tip, with PC connector

Vented tip



<sup>\* 5</sup>pcs/box, 100pcs/ctn

<sup>\*</sup> Luer-lcok & 15mm PC Connector

We have many ranges of Endotracheal tube Bougies.

# 7. Tracheal Tube Introducer (Bougies), Vented tip, with 2 connectors


Vented tip



<sup>\* 5</sup>pcs/box, 100pcs/ctn

## 8. Tracheal Tube Introducer (Exchange Tube)

Vented tip



#### **FEATURES AND BENEFIT:**

- Used to facilitate uncomplicated, atraumatic endotracheal tube exchange.
- Made of medical grade PE, LATEX Free, Flexible.
- Blunt tip is atraumatic to internal structures.
- Distal side port help to ensure sufficient air flow.
- 1 CM markings facilitate accurate placement with endotracheal tubes.
- Removable adapter permit ventilatory if necessary during exchange procedure.

REF number	Size	OD	Length	Recommended ETT size
PA2107D08	8Fr	2.7mm	53.5cm	ID 3.0mm - 4.0mm
PA2107D11	11Fr	3.7mm	80cm	ID 4.5mm - 6.0mm
PA2107D14	14Fr	4.7mm	80cm	ID 6.0mm - 10.0mm
PA2107D19	19Fr	6.3mm	80cm	ID 8.0mm - 10.0mm

<sup>\* 5</sup>pcs/box, 100pcs/ctn



\* Luer-lock & 15mm PC Connector



## 9. Connectors and Adapters

\* Our removable Adapters can be attached to Bougies and connect with ventilation devices.



15mm TPE connector

REF number	Material	Suitable Bougie size
PA2103TPE1	TPE	6Fr, 8Fr, 10Fr, 11Fr
PA2103TPE2	TPE	14Fr, 15Fr



Luer-lock & 15mm PC connector

REF number	Material	Suitable Bougie size
PA2103PC1	PC	8Fr, 10Fr, 11Fr
PA2103PC2	PC	14Fr, 15Fr, 19Fr

#### **ABOUT US:**

FORMED is a professional factory focus on medical devices located in Hangzhou, Zhejiang province, China. We have more than 20 years experience in Medical Devices, including design, development, regulation, production, quality control, sales and services.


Our main products cover Anaesthesia, Respiratory, Gastroenterology, Urology, Wound Drainage and Surgical. Like: Laryngeal Masks, Introducers&Bougies, Gastric Calibration Tubes, Tracheal Tubes, Bronchial Tubes, Oxygen Masks, Resuscitators, CPAP Masks, Stylets, Circuits&Filters, Suction, Catheters, Enteral Feeding etc. We also cooperate with customers to design new products.

We have ISO 13485, MDR CE, US FDA, UK MHRA, Australia TGA, Brazil ANVISA, Saudi Arabia SFDA certificates. We always keep providing high quality products to customers.



































### HANGZHOU FORMED MEDICAL DEVICES CO.,LTD.

Add: No.292 Gonghe North Road, Xindeng Town, Fuyang District, Hangzhou, Zhejiang Province, China, 311404 Email: info@formedtech.com

#### 杭州福美医疗器械有限公司

浙江省杭州市富阳区新登镇共和北路292号

#### **Sales Center**

Add: No. 3-014, Zhongnan Garden Business Building, Wuyi County, Jinhua City, Zhejiang Province, China, 321200 Tel: +86-579-8787 6372 Ext 8008



WhatsApp

www. formed tech. com

Sky Gao | Sales Dept. Manager 高 俊 | 销售部经理

Mobile: +86-183 9599 0577 Email: sky@formedtech.com