



# MC-150 Cryoset LN<sub>2</sub> Unit

---




**CRYOMED BIOMEDICAL**



## Spray tip\_A

**Code no.: 150.101**


Gold plated spray tip-1mm



## Spray tip\_B

**Code no.: 150.102**


Gold plated spray tip-0,8mm



## Spray tip\_C

**Code no.: 150.103**


Gold plated spray tip-0,6mm



## Spray tip\_D

**Code no.: 150.104**


Gold plated spray tip-0,4mm



## Spray tip\_E

**Code no.: 150.105**

Gold plated spray tip-0,3mm




## Conical Spray tip\_A

---

**Code no.:** 150.106

Gold plated spray tip-1mm




## Conical Spray tip\_B

---

**Code no.:** 150.107

Gold plated spray tip-0,8mm




## Conical Spray tip\_C

---

**Code no.:** 150.108

Gold plated spray tip-0,6mm




## Conical Spray tip\_D

---

**Code no.:** 150.109

Gold plated spray tip-0,4mm




## Acne Spray tip

---

**Code no.:** 150.110


Gold plated spray tip  
Maximizes the effect in acne interventions



## Cryo chamber\_6

**Code no.: 150.111**

Gold plated cryo chamber-6mm  
It is designed to perform a pool effect on deep lesions



## Cryo chamber\_10

**Code no.: 150.112**

Gold plated cryo chamber-10mm  
It is designed to perform a pool effect on deep lesions



## Cryo chamber\_15

**Code no.: 150.113**

Gold plated cryo chamber-15mm  
It is designed to perform a pool effect on deep lesions



## Cryo chamber\_18

**Code no.: 150.114**

Gold plated cryo chamber-18mm  
It is designed to perform a pool effect on deep lesions



## Cryo chamber\_25

**Code no.: 150.115**

Gold plated cryo chamber-25mm  
It is designed to perform a pool effect on deep lesions



## Cryo chamber\_30

---

**Code no.: 150.116**

Gold plated cryo chamber-30mm  
It is designed to perform a pool effect on  
deep lesions



## Contact probe tip\_0,5

---

**Code no.: 150.117**

Gold plated tip-0,5mm  
In mini lesions, treatment is carried out  
by contact



## Contact probe tip\_1

---

**Code no.: 150.118**

Gold plated tip-1mm  
In mini lesions, treatment is carried out  
by contact




## Contact probe tip\_2

---

**Code no.: 150.119**

Gold plated tip-2mm  
In mini lesions, treatment is carried out  
by contact




## Contact probe tip\_3

---

**Code no.: 150.120**

Gold plated tip-3mm  
In mini lesions, treatment is carried out  
by contact



### Contact probe tip\_4

---

**Code no.: 150.121**

Gold plated tip-4mm  
In mini lesions, treatment is carried out  
by contact




### Contact probe tip\_5

---

**Code no.: 150.122**

Gold plated tip-5mm  
In mini lesions, treatment is carried out  
by contact



### Contact probe tip\_6

---

**Code no.: 150.123**

Gold plated tip-6mm  
In mini lesions, treatment is carried out  
by contact




### Contact probe tip\_8

---

**Code no.: 150.124**

Gold plated tip-8mm  
In mini lesions, treatment is carried out  
by contact




### Contact probe tip\_10

---

**Code no.: 150.125**

Gold plated tip-10mm  
In large lesions, treatment is carried out  
by contact



## Contact probe tip\_12

---

**Code no.: 150.126**

Gold plated tip-12mm  
In large lesions, treatment is carried out  
by contact



## Contact probe tip\_15

---

**Code no.: 150.127**

Gold plated tip-15mm  
In large lesions, treatment is carried out  
by contact



## Contact probe tip\_20

---

**Code no.: 150.128**

Gold plated tip-20mm  
In large lesions, treatment is carried out  
by contact




## Contact probe tip\_25

---

**Code no.: 150.129**

Gold plated tip-25mm  
In large lesions, treatment is carried out  
by contact




## Long-Bent Spray tip

---

**Code no.: 150.130**

Gold plated spray tip  
For lesions that are difficult to reach;  
14-16-18-20G sizes




## Long Spray tip

---

**Code no.: 150.131**

Gold plated spray tip  
For lesions that are difficult to reach;  
14-16-18-20G sizes




## Long Spray tip\_2

---

**Code no.: 150.132**

Gold plated spray tip-2cm  
For lesions that are difficult to reach;  
14-16-18-20G sizes




## Long Spray tip\_3

---

**Code no.: 150.133**

Gold plated spray tip-3cm  
For lesions that are difficult to reach;  
14-16-18-20G sizes




## Long Spray set

---

**Code no.: 150.134**

Gold plated spray tips-4pcs  
For lesions that are difficult to reach;  
14-16-18-20G sizes




## Angle Extension adapter\_90

---

**Code no.: 150.135**

Gold plated angle extension adapter-90°  
It allows the device to stay in an upright position by changing the spray angle during the cryotherapy intervention.  
Suitable for use with any Bimed® spray/contact tips and cryo chambers.






## Angle Extension adapter\_45

**Code no.: 150.136**


Gold plated angle extension adapter-45°  
It allows the device to stay in an upright position by changing the spray angle during the cryotherapy intervention.  
Suitable for use with any Bimed® spray/contact tips and cryo chambers.



## Flexible Extension adapter

**Code no.: 150.137**


Gold plated flexible extension adapter  
It allows the device to stay in an upright position by changing the spray angle during the cryotherapy intervention.  
Suitable for use with any Bimed® spray/contact tips and cryo chambers.



## Extension adapter

**Code no.: 150.138**


Gold plated extension adapter  
It allows the device to stay in an upright position by changing the spray angle during the cryotherapy intervention.  
Suitable for use with any Bimed® spray/contact tips and cryo chambers.



## Extension adapters set

**Code no.: 150.139**


Gold plated extension adapters-4pcs  
It allows the device to stay in an upright position by changing the spray angle during the cryotherapy intervention.  
Suitable for use with any Bimed® spray/contact tips and cryo chambers.



## Back Vent adapter

**Code no.: 150.140**

Gold plated back vent adapter  
Prevents intermittent spraying and clogging. Increases liquid nitrogen flowing.



## Luer-Lock adapter

---

**Code no.: 150.141**

**Gold plated Luer-Lock adapter**  
Used to spray the liquid nitrogen with the disposable needles.




## Cryo Cones

---

**Code no.: 150.142**

**Cryo cones**  
Used in irregularly shaped lesions to take the shape of the lesion, to limit the spray area and to concentrate the liquid nitrogen.




## Cryo Plate\_1,5

---

**Code no.: 150.143**

**Transparent Cryo plate-1,5mm**  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas.




## Cryo Plate\_3

---

**Code no.: 150.144**

**Transparent Cryo plate-3mm**  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas.




## Cryo Plate\_5

---

**Code no.: 150.145**

**Transparent Cryo plate-5mm**  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas.




## Cryo Plate\_8

---

**Code no.: 150.146**

Transparent Cryo plate-8mm  
Collects the spray in a certain area by  
limiting the lesion in lesions on the  
sensitive areas.




## Cryo Plate\_10

---

**Code no.: 150.147**

Transparent Cryo plate-10mm  
Collects the spray in a certain area by  
limiting the lesion in lesions on the  
sensitive areas.




## Cryo Plate\_12

---

**Code no.: 150.148**

Transparent Cryo plate-12mm  
Collects the spray in a certain area by  
limiting the lesion in lesions on the  
sensitive areas.




## Cryo Plate\_15

---

**Code no.: 150.149**

Transparent Cryo plate-15mm  
Collects the spray in a certain area by  
limiting the lesion in lesions on the  
sensitive areas.




## Cryo Plate\_18

---

**Code no.: 150.150**


Transparent Cryo plate-18mm  
Collects the spray in a certain area by  
limiting the lesion in lesions on the  
sensitive areas.



## Mini Cryo Plate set

**Code no.: 150.151**


Transparent Cryo plates-4pcs  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas.



## Maxi Cryo Plate set

**Code no.: 150.152**


Transparent Cryo plates-4pcs  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas.



## Cryo Plate set

**Code no.: 150.153**

Transparent Cryo plates-8pcs  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas.



## Cervical probe\_0519

**Code no.: 150.154**

Gold plated gynecological probe-0519



## Cervical probe\_1905

**Code no.: 150.155**

Gold plated gynecological probe-1905



## Cervical probe\_1910

**Code no.: 150.156**

Gold plated gynecological probe-1910



## Cervical probe\_1920

**Code no.: 150.157**

Gold plated gynecological probe-1920



## Cervical probe\_2500

**Code no.: 150.158**

Gold plated gynecological probe-2500



## Cervical probe\_2507

**Code no.: 150.159**

Gold plated gynecological probe-2507



## Protective goggles

**Code no.: 150.160**

Protective goggles  
Used to protect the eyes during liquid nitrogen  
filling process.



## Protective gloves

---

**Code no.: 150.161**

**Protective gloves**  
Used to protect the hands during liquid nitrogen filling process.



## Cryo Plate\_4 holes

---

**Code no.: 150.162**

**Transparent Cryo plate with 4 holes**  
Collects the spray in a certain area by limiting the lesion in lesions on the sensitive areas. According to the lesion size, holes of four different diameters are preferred.

