TECHARTMED ®

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery

ET500T Electric Operating Table

Developments made in the conjunction with doctors: this is one of the basic principles which has made TECHARTMED the world market leader for operating tables. Close cooperation and years of experience have resulted in tried and tested products which facilitate the work in the operating room and help to increase the safety and comfort of the

Best example: this table can withstand even extreme loads in any position. The structure enables it to be adapted to match the size of the patient and in no time at all it can be equipped for any surgical discipline. A longitudinal shift of up 350mm makes the operating table even more versatile.

Your Ideal Operating Table Supplier



NINGBO TECHART MEDICAL EQUIPMENT CO.,LTD

Add: No.99 South Hongtang Road Ningbo China 315033

Tel: +86-574-83076915

Fax: +86-574-87137656

E-mail: director@techartmed.com sales01@techartmed.com sales02@techartmed.com sales03@techartmed.com

Web: www.techartmed.com



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery



NINGBO TECHART MEDICAL EQUIPMENT CO., LTD

TECHARTMED®

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery

Electric Operating Table



Engineered for increasingly complex surgical scenes, ET500T comes with intelligent automatic control system satisfying various demands of general surgeries. Extra wide tabletop, extreme loading system, long horizontal sliding that can be suitable for both X-ray and C-arm. ET500T is an ideal electric operating table.

- Longitudinal slide up to 350mm for free access to C arm
- Extreme weight load capacity of 250kgs with highest safety and stability
- Memory pad with antistatic ,waterproof design
- Standard override panel control makes the control safer
- Battery power for 50–80 operations

PARAMETER ET500T Electric Operating Table

Move	Model:ET500T			
Tableton trandelanhurs	Reverse Trendelenburg	≥25°		
Tabletop trendelenburg	Trendelenburg	≥25°		
Lateral Tilt (L	eft And Right)	≥15°		
	Up	≥40°		
Head plate	Down	≥90°		
	Flexible	≥40°		
Back plate	Up	≥80°		
	Down	≥20°		
Leg plate	Up	≥15°		
	Down	≥90°		
	Outward	≥90°		
Kidney B	ridge Up	110mm		
Horizont	≥350mm			



1, Longitudinal shift: the longitudinal shift of up to 350MM enables optimum access with the C-arm without needing to reposition the patient or to reverse the positioning of the patient. Each selected fluoroscopy point may easily be moved in front of or behind the column and ensures comfortable access to the patient from initiation to

- 2, It can easily be moved to any desired position as four castors.
- 3, Leg plate: imported gas spring for easy controls.
- 4, Two joint head plate.
- 5, Built-in kidney bridge.
- 6, Double control system.





Well Sold For 40 Years Worldwide Meet The Basic Demands Of General Surgery

Your Ideal Operating Table Supplier

Model Site	Length	Width	Tabletop lowest height	Elevation range
ET500T	≥2100mm	550mm	≤700mm	300mm

Anaesthesia screen	1piece	Body support	1 pair	Leg support	1 pair	Hand remote	1 piece
Kidney bridge handle	1 piece	Arm support	1 pair	Foot support	1 pair	Mattress	1 set
Shoulder support	1 pair	Long fixing clamp	2 pieces	Fixing clamp	7 pieces	Power Line	1piece





