

PT Oneject Indonesia is the only Smart Syringe manufacturer in Indonesia with a certificate of Domestic Component Level (TKDN or Tingkat Komponen Dalam Negeri) approximately 60% from Indonesia's Ministry of Industry.

PT Oneject Indonesia also has complete ISO Certificates (ISO 9001:2015, ISO 13485:2016), CE Certificate from TUV Nord Indonesia, GMP, and Prequalification PQS-WHO. In the future, PT Oneject Indonesia will always promote local products in the global industry.



One Injection, One Syringe

PT Oneject Indonesia Sentul Factory:

Jl. Olympic Raya No.B9, Kawasan Industri
Sentul, Bogor 16810, Indonesia

Bekasi Factory:

Jl. Anggrek VII No.41, Kawasan Industri Terpadu
Indonesia China (KITIC), Kab. Bekasi 17337, Indonesia

Tel. : +62 21 8791 7422 (Hunting)

Fax. : +62 21 8790 2761

E-mail : marketing@oneject.co.id



 **oneject**

Auto Disable Syringe

Every man,
woman & child
deserves
the right to a
safe injection

www.oneject.co.id

Oneject Auto Disable Syringe

ONEJECT Auto Disable Syringe design uses a simple yet effective technology by Star Syringe Limited (UK). After injection, the plunger will lock and immobilize automatically making the syringe unusable. This design is WHO / PQS approved and certified to all Auto Disable Syringes from PT Oneject Indonesia.



AUTO DISABLE SYRINGE ADVANTAGES:

1. Auto Disable Syringe can only be used once.
2. After injection, plunger automatically locked and will break when retracted.
3. Minimizes disease contamination from syringes.
4. Minimizes syringes misuse such as use more than once.
5. Protect Patients and Healthcare Workers from contamination .
6. Increase patient trust in Healthcare Workers.



Oneject AD Syringe 0.05ml
26G x 3/8" (10 mm)
PQS-WHO No: E008/030



Oneject AD Syringe 0.3ml
24G x 1" (25 mm)
PQS-WHO No: E008/086



Oneject AD Syringe 0.5ml
24G x 1" (25 mm)
PQS-WHO No: E008/019



Oneject Hypodermic Syringe with Reuse Prevention Feature 1 ml
23G x 1" (25 mm)
PQS-WHO No: E013/023



Oneject Hypodermic Syringe with Reuse Prevention Feature 2 ml
23G x 1 1/4" (32 mm)
PQS-WHO No:



Oneject Hypodermic Syringe with Reuse Prevention Feature 3 ml
23G x 1 1/4" (32 mm)
PQS-WHO No: E013/019



Oneject Hypodermic Syringe with Reuse Prevention Feature 5 ml
21G x 1 1/2" (38 mm)
PQS-WHO No: E013/024



**one injection,
one syringe**