
## FLAT STERILIZATION ROLLS

- Sudemed sterilization rolls are produced from transparent multi-layer PET/PP film and medical kraft paper.
- Optional sterilization rolls can be produced in green, blue, pink, purple, orange or any other color you request.
- It is manufactured from the world's highest quality 60 gr. and 70 gr. medical kraft papers.
- Sudemed sterilization rolls are produced without the use of solvents in high-tech production processes that prevent the risk of tearing.
- All rolls are available with water-based, non-toxic vapor and ethylene oxide (EO) process indicators. An optional formaldehyde indicator is also available.
- Conforms to EN 868-5 standard



Product Code	Product Name	Size	Pack
SDM05200	Sterilization Roll	5 cm x 200 m	12
SDM75200	Sterilization Roll	7,5 cm x 200 m	8
SDM10200	Sterilization Roll	10 cm x 200 m	8
SDM15200	Sterilization Roll	15 cm x 200 m	4
SDM20200	Sterilization Roll	20 cm x 200 m	4
SDM25200	Sterilization Roll	25 cm x 200 m	2
SDM30200	Sterilization Roll	30 cm x 200 m	2
SDM35200	Sterilization Roll	35 cm x 200 m	2
SDM40200	Sterilization Roll	40 cm x 200 m	2
SDM50200	Sterilization Roll	50 cm x 200 m	2
SDM60200	Sterilization Roll	60 cm x 200 m	2







## STERILIZATION POUCH

- Sudemed sterilization pouches, ready for packaging, closed on three sides, are produced from transparent, multi-layer PET/PP film and medical kraft paper.
- Optional sterilization pouches can be produced in green, blue, pink, purple, orange or any other color you request.
- Manufactured from the world's highest quality 60 gr and 70 gr medical kraft paper.
- In the latest technological production processes, it consists of reinforced film that prevents the risk of tearing. Solvent is not used in our film laminations.
- All pouches have water-based, non-toxic vapor and ethylene oxide (EO) process indicators. Formaldehyde indicators are also available upon request.
- Conforms to EN 868-5 standard.



## TYVEK® STERILIZATION POUCHES

- Ready for packaging, Sudemed sterilization pouches made of three-side sealed Tyvek® material are made of transparent, multi-layer PET/PE film and Tyvek® material with high tear resistance
- In addition to being used in hydrogen peroxide gas plasma sterilization in hospitals, it is also used in EO and Gamma sterilization in the industrial field.
- It consists of a reinforced film that prevents the risk of tearing in the latest technological production processes. Solvents are not used in our film laminations.
- All pouches contain water-based, non-toxic hydrogen peroxide (H2O2) process indicators.
- Complies with EN 868-5 standard.