



Patient Monitor

Multi-use smart monitor's 7-inch model

IP-3010



7-inch High
resolution screen





ECG(3 or 5 lead), SpO2, NIBP,
2-IBP, 2 temperature,
respiration (Opt.EtCO2/printer)



Perfusion range
0.05%~20%



Advanced SpO2 technology
(Motion artifact)





7-inch High resolution screen

The system displays 7 parameters and 5 waveforms on a 7-inch wide screen (800×480). Each waveform is presented in high-resolution images with clear, identifiable colors, enhancing user comfort with TFT color and an intuitive user interface.


Thermal printer

It comes with a built-in thermal printer, allowing for immediate data printing without the need for additional connections or devices.




ECG analysis

It detects premature ventricular contractions and 12 types of arrhythmias through ECG recording and waveform interval measuring.



ETCO2 Module

CO₂ pressure is measured directly by inserting sensors into the intubation site or by sidestream measurement. The EtCO₂ Module (Sidestream) can be optionally integrated and used without additional connections or attachments.



Nurse call

It quickly informs medical staff when a patient's abnormal symptoms are detected, allowing for a rapid response in case of an emergency.

I Manufacturer



Infunix Technology Co., Ltd.

T +82-2-3448-2800

F +82-2-546-8346

E sales@infunix.com

I Sales





IP-4050

PureScope series

Bedside monitor



7 Parameter
6 waveforms



10.4 inch
touch screen



ECG (3 or 5 lead), SpO2, NIBP,
2-IBP, 2-temperature, respiration
(Opt.EtCO2 / printer)




Detects 12 types of
ECG arrhythmia and
displays PVC and ST level



Notification system when
it detects an abnormality
in a patient



Central monitoring system
(Monitoring 32 beds)



ECG analysis

- 12 Kinds of arrhythmia
- Arrhythmia alarm
- Display of PVC & ST level
- Pacemaker detection

Trend

72 hours of graphic and tabular trend

Built-in EtCO2 module

- Easy connect EtCO2 cannula to the monitor directly
- Convenient & user-friendly use
- Prevent from destroy or loss of external module

Optional

- Touch screen
- Rolling stand & wall mount
- Printer



Central monitoring system

Bidirectional communication with patient monitors
Monitoring 4, 8, 16, 32 beds
Waveform, parameter and alarm are saved.
Support hl7

| Manufacturer



Infunix Technology Co., LTD

T +82-2-3448-2800
F +82-2-546-8346
E sales@infunix.com

| Sales



GENORAY Co.,Ltd.

T +82-31-5178-5500
E genoray@genoray.com
W www.genoray.com

GENORAY America Inc.

T +1-855-436-6729
E inquiry@genorayamerica.com
W www.genorayamerica.com

GENORAY EU GmbH

T +49-30-50969498
E info@genorayeuropa.com
W www.europe.genoray.com

GENORAY Japan

T +81-45-620-4971
E Info@genorayjapan.co.jp
W www.genorayjapan.co.jp

GENORAY Türkiye

T +90-216-463-57-37
+90-216-463-57-53
E genorayturkiye@genoray.com


GENORAY Shanghai

T +86-021-6785-8538
W www.genoray.com.cn

IP-5000

PureScope series

Multi parameter patient monitor



7 Parameter
9 waveforms



15 Inch
touch screen



Configurable color of
parameters & waveforms





Detects 12 types of
ECG arrhythmia and
displays PVC and st level



Operates 4 hours
charged battery



Central monitoring system
(Monitoring 32 beds)



ECG analysis

- 12 Kinds of arrhythmia
- Arrhythmia alarm
- Display of PVC & ST level
- Pacemaker detection

Trend

168 Hours of graphic and tabular trend

SpO2 averaging (4/8/12/16)

Alarm duration adjustable

Optional

- Touch screen
- Rolling stand & wall mount
- Printer



Central monitoring system

Bidirectional communication with patient monitors
Monitoring 4, 8, 16, 32 beds
Waveform, parameter and alarm are saved.
Support hl7

| Manufacturer



Infunix Technology Co., LTD

T +82-2-3448-2800
F +82-2-546-8346
E sales@infunix.com

| Sales



GENORAY Co.,Ltd.

T +82-31-5178-5500
E genoray@genoray.com
W www.genoray.com

GENORAY America Inc.

T +1-855-436-6729
E inquiry@genorayamerica.com
W www.genorayamerica.com

GENORAY EU GmbH

T +49-30-50969498
E info@genorayeuropa.com
W www.europe.genoray.com

GENORAY Japan

T +81-45-620-4971
E Info@genorayjapan.co.jp
W www.genorayjapan.co.jp

GENORAY Türkiye

T +90-216-463-57-37
+90-216-463-57-53
E genorayturkiye@genoray.com

GENORAY Shanghai

T +86-021-6785-8538
W www.genoray.com.cn