

PRODUCT LIST



AGM-513S

GlucoDr.S



AGM-4000

GlucoDr.auto™



AGM-2200

GlucoDr
SUPER SENSOR



AGM-2100

GlucoDr™

As a promising research-focused company, **Allmedicus** has continued to develop and supply innovative diagnostic medical instruments since 1997. **Allmedicus** aims to improve the quality of human healthcare by developing accurate, efficient, technologically advanced and cost-effective monitoring system based on biotechnology and medical engineering for both individuals and healthcare professionals.

GlucoDr.S

Blood Glucose Monitoring System

Specifications

Product Name	GlucoDr.S™ Blood Glucose Monitoring System
Model Name	AGM-513S
Assay Method	Electrochemical method
Sample Type	Fresh capillary and venous whole blood
Sample Volume	0.5 µl
Measurement Range	20~600 mg/dL (1.1~33.3 mmol/L)
Measurement Time	5 seconds
Calibration	Plasma-equivalent
Battery Type	One 3-volt lithium battery (coin cell type CR2032)
Battery Life	Approximately 1,000 tests
Unit of Measurement	mg/dL or mmol/L
Operating Ranges	Temperature: 10~40°C (50~104°F) Relative Humidity: 15~85% Altitude: Up to 3048m (10,000 feet) Hematocrit: 20~65%
Dimension	50.1 X 87.1 X 17.9 (mm), 2.0 X 3.4 X 0.70 (inches)
Weight	Approximately 40g or 1.4 oz. (with battery)
Display	44.0 X 38.0 (mm) LCD, 1.7 X 1.5 (inches) LCD
Memory Capacity	500 test results (with date, time & event)
PC Interface	Micro USB
Automatic Shutoff	5 minutes after last action

All Medicus is always nearby you.

 **allmedicus**
www.allmedicus.com

All Medicus Makes
Standard for Healthcare


GlucoDr.S

Blood Glucose Monitoring System




All New ISO Compliance!!
ISO 15197:2013 (EN ISO 15197:2015)

AGM-513S



 **allmedicus**
www.allmedicus.com

AGM-513S TEST STRIP



Glucose Dehydrogenase (GDH-FAD)


Gold Electrode



GDH-FAD

Superior accuracy
Interference-free from non-glucose sugar


AGM-513S TEST STRIP




Super sip-in technology

Fast and easy blood absorption

Gold Electrode



Absorbent Hole



Side blood application

Avoid touching enzyme chamber with sweaty fingers to assure accurate testing

High quality & Handy BGMS

AGM-513S




GlucoDr.S


Blood Glucose Monitoring System



- No Coding
- 0.5 ul / 5 sec
- Glucose Dehydrogenase (GDH-FAD)
- Wide Hematocrit Range (20%~65%)
- Strip Release Button
- High & Low Level Indicator



No Coding




Release Button



NFC / Bluetooth (Option)



99.99% Gold Electrode



Micro USB

Manufactured by
All Medicus Co., Ltd.

#7102~7107, #7402, #7403, #7406, 140, Beolmal-ro, Dongan-gu,
Anyang-si, Gyeonggi-do 14057 REPUBLIC OF KOREA
www.allmedicus.com