

THE TOP IN QUALITY for all areas of healthcare

TOP PLASTIC SYRINGE

An easy-to-use syringe which is highly airtight and smooth-working



Capacity	Without Needle		With Needle		Package	
1mL	4			25g × 25mm	Without Needle : 100 pcs./box 24 boxes/carton With Needle : 100 pcs./box 12 boxes/carton	
1mL (for Tuberculin)				27G × 13 mm	100 pcs./box 18 boxes/carton	
2.5mL		٩///		23G × 25mm	100 pcs./box 12 boxes/carton	
5mL		٩///		22g × 32mm	100 pcs./box 12 boxes/carton	
10mL		٩///		21g × 38mm	100 pcs./box 8 boxes/carton	
20mL		٩///		$21 \text{G} \times 38 \text{mm}$	50 pcs./box 8 boxes/carton	
30mL		٩///			50 pcs./box 8 boxes/carton	
50mL		٩///	25 pcs./box 12 boxes/carto		25 pcs./box 12 boxes/carton	

TOP Corporation

19-10 Senjunakai-cho, Adachi-ku, Tokyo 120-0035, Japan Website : http://top-tokyo.co.jp



TOP INJECTION NEEDLE

Sharp needle points provide smooth injections and lessen patient anxiety



The needle hub is highly transparent, making it easy to check flashback when puncturing a blood vessel. Each size is color-coded according to the ISO 6009, and the size can be identified at a glance.

Gauge / Color	Needle length	Package	Gauge / Color	Needle length	Package
18G (1.27mm)	11/2(38mm)		23G	1 (25mm)	100 pcs. / Box 30 Box / Carton
	1/2(13mm)		(0.63mm)	5/8(16mm)	
19G (1.10mm)	11/2(38mm)		24G	1 (25mm)	
	5/8(16mm)	100 peo / Dev	(0.55mm)	5/8(16mm)	
21G (0.80mm)	5/8(16mm)	100 pcs. / Box 30 Box / Carton	25G(0.50mm) -	5/8(16mm)	
	, , , ,		250(0.501111)	1(25mm)	
	1/2(13mm)		26G (0.45mm)	1/2(13mm)	
22G (0.70mm)	11/4(32mm)		27G	3/4(19mm)	-
	5/8(16mm)		(0.40mm)	1/2(13mm)	

TOP Corporation

19-10 Senjunakai-cho, Adachi-ku, Tokyo 120-0035, Japan Website : https://www.top-tokyo.co.jp